

## **SUPPLEMENTAL TECHNICAL BULLETIN**

### **ST – 2006 – 18**

**Title:** Importing PAP-8E Aggregation Data into Excel

**Bulletin No:** ST-2006-18


This Supplemental Technical Bulletin (ST) has been developed as a laboratory aid. This ST does not alter, revise or change the information provided in the Technical Bulletin included with each product. In accordance with Good Laboratory Practice and regulatory requirements, each laboratory must develop, validate and adopt its own written procedures.

---


### Instructions for Importing PAP-8E Aggregation Data into Excel

1. Open the supplied Excel Aggregation Graph Template 1 attached with this bulletin. Go to File, then Open
2. Look in: Select Local Disk (C)
3. Select PAP-8E\_Export, then open
4. Files of type, Select All Files
5. File name, Click on the selected exported file, then open
6. Follow screen instructions to import file to this Excel Worksheet  
Text Import Wizard  
Step 1 – Check “delimited”, Next  
Step 2 – Check delimiters, “Tab, Semicolon, Comma” – Next  
Step 3 – Finish
7. Find, Select and Copy (Control C) the data that you want to graph
8. Paste (Control V) the selected data into the (Excel Aggregation Graph Template 1) worksheet. Place the cursor in Cell 1A, then paste (Control V).
9. Keep a clean copy of the “Excel Aggregation Graph Template 1” for each use.

vW Ristocetin CoFactor Assay, 4 point curve, data imported into Excel:



	A	B	C	D	E	F	G	H	I	J	K	L
1	3.762	6.818	2.888	-0.812								
2	4.369	5.195	3.04	-1.299								
3	3.762	2.002	2.026	-1.299								
4	3.236	-0.487	4.407	-1.677								
5	4.571	-0.325	7.548	0.541								
6	7.241	1.515	11.803	3.084								
7	9.628	2.597	11.905	5.519								
8	9.223	2.327	11.854	5.898								
9	9.79	2.597	10.132	6.115								
10	9.426	4.6	9.98	5.682								
11	10.639	7.359	5.941	2.246								
12	6.388	8.929	6.273	-0.952								
13	3.683	8.117	3.863	-2.6								
14	-1.011	8.171	4.525	-1.537								
15	-1.149	5.47	2.483	-0.355								
16	-1.425	4.06	3.366	-0.946								
17	-1.287	1.705	3.091	-0.827								
18	-1.241	2.175	3.863	-3.546								
19	-2.483	1.529	3.863	-3.073								
20	-2.897	1.529	4.249	-2.66								
21	-3.448	0.823	4.857	-2.364								
22	-2.345	0.647	4.249	-2.009								
23	-2.667	-0.235	3.918	-3.487								



	A	B	C	D	E	F	G	H	I	J	K	L
1												
2	Seconds	Minutes		Channe 1	Channel 2	Channel 3	Channel 4	Channel 5	Channel 6	Channel 7	Channel 8	
3	0.0	0.000		3.762	6.818	2.888	-0.812	0	0	0	0	
4	0.5	0.008		4.369	5.195	3.04	-1.299	0	0	0	0	
5	1.0	0.017		3.762	2.002	2.026	-1.299	0	0	0	0	
6	1.5	0.025		3.236	-0.487	4.407	-1.677	0	0	0	0	
7	2.0	0.033		4.571	-0.325	7.548	0.541	0	0	0	0	
8	2.5	0.042		7.241	1.515	11.803	3.084	0	0	0	0	
9	3.0	0.050		9.628	2.597	11.905	5.519	0	0	0	0	
10	3.5	0.058		9.223	2.327	11.854	5.898	0	0	0	0	
11	4.0	0.067		9.79	2.597	10.132	6.115	0	0	0	0	
12	4.5	0.075		9.426	4.6	9.98	5.682	0	0	0	0	
13	5.0	0.083		10.639	7.359	5.941	2.246	0	0	0	0	
14	5.5	0.092		6.388	8.929	6.273	-0.952	0	0	0	0	
15	6.0	0.100		3.683	8.117	3.863	-2.6	0	0	0	0	
16	6.5	0.108		-1.011	8.171	4.525	-1.537	0	0	0	0	
17	7.0	0.117		-1.149	5.47	2.483	-0.355	0	0	0	0	
18	7.5	0.125		-1.425	4.06	3.366	-0.946	0	0	0	0	
19	8.0	0.133		-1.287	1.705	3.091	-0.827	0	0	0	0	
20	8.5	0.142		-1.241	2.175	3.863	-3.546	0	0	0	0	
21	9.0	0.150		-2.483	1.529	3.863	-3.073	0	0	0	0	
22	9.5	0.158		-2.897	1.529	4.249	-2.66	0	0	0	0	
23	10.0	0.167		-3.448	0.823	4.857	-2.364	0	0	0	0	

