

SUPPLEMENTAL TECHNICAL BULLETIN ST – 2011 – 02

Title: Summary of ISTH 2009 Meeting - Updated Guidelines for Lupus Anticoagulant Detection

Bulletin No: ST-2011-02

This Supplemental Technical Bulletin (ST) has been developed as a laboratory aid. This ST does not alter, revise or change the information provided in the Technical Bulletin included with each product. In accordance with Good Laboratory Practice and regulatory requirements, each laboratory must develop, validate and adopt its own written procedures.

Summary of ISTH 2009 Meeting - Updated Guidelines for Lupus Anticoagulant Detection

At the ISTH meeting in Boston on July 11, 2009, the Scientific and Standardization Committee on Lupus Anticoagulant, chaired by Vittorio Pengo, MD, of the University Of Padua School of Medicine, Dr. Pengo and distinguished members of the committee presented an update of the Lupus Anticoagulant (LA) Detection Guidelines, first published in 1995.

Here are some of the highlights.

- Appropriateness for testing is graded low: arterial or venous thromboembolism (VTE) in elderly patients; moderate: unexpected prolonged partial thromboplastin time (PTT) in asymptomatic subjects, recurrent spontaneous early pregnancy loss or provoked VTE in patients <50 years of age; and high: unprovoked VTE or arterial thrombosis <50 years of age, thrombosis at unusual sites, late pregnancy loss, any thrombosis or pregnancy morbidity in patients with autoimmune diseases.
- Double centrifuge to ensure platelet poor plasma (PPP, <10,000/uL), quick-freeze if testing is delayed, thaw at 37°C.
- Test using two clot-based assays employing separate clotting principles. The dilute Russell Viper Venom Time assay is the most robust. The other may be low-phospholipid PTT with silica activator. PTTs with kaolin activator are problematic in automated coagulometers, ellagic acid is insensitive to LA and both are contraindicated.
- The following tests are not recommended: Dilute Prothrombin Time, Ecarin and Textarin times, and Kaolin Clotting Time.

- Perform a thrombin time to detect therapeutic heparin and absorb heparin or collect at another time.
- For mixing studies, pooled normal plasma (PNP) is prepared “ad hoc” (home-made) by double centrifugation to ensure it is PPP. It must provide 100% activity for all clotting factors. Commercial lyophilized or frozen normal plasmas can be used if they fulfill these specifications.
- For mixing studies, the cut-off is the 99th percentile of the distribution or an Index of Circulating Anticoagulant (ICA): computed as follows...

$$\text{ICA}\% = [(\text{clot time of patient-PNP mixture} - \text{clot time of PNP})/\text{clot time of patient}] \times 10$$

- The confirmatory test may use bilayer or hexagonal phase phospholipid. Freeze/thawed platelets are not recommended because of poor batch-to-batch consistency.
- Confirmatory test cut-off is defined as follows...

$$[(\text{screen} - \text{confirm})/\text{screen}] \times 100$$

- Cut-off values are developed using the local reagent/coagulometer combination on at least 40 adult healthy donors <50 years of age.

References

1. Pengo V, Tripodi A, Reber G, Rand JH, Ortel TL, Galli M, de Groot PG. Update of the guidelines for lupus anticoagulant detection. *J Thromb Haemost* 2009; 7: 1737–40.
2. Brandt JT, Triplett DA, Alving B, Scharrer I. Criteria for the diagnosis of lupus anticoagulants: an update. On behalf of the Subcommittee on Lupus Anticoagulant/Antiphospholipid Antibody of the Scientific and Standardisation Committee of the ISTH. *Thromb Haemost* 1995; 74:1185–90.
3. Miyakis S, Lockshin MD, Atsumi T, Branch DW, Brey RL, Cervera R, Derksen RHWM, de Groot PG, Koike T, Meroni PL, Reber G, Shoenfeld Y, Tincani A, Vlachoyiannopoulos PG, Kritis SA. International consensus statement on an update of the classification criteria for definite antiphospholipid syndrome (APS). *J Thromb Haemost* 2006; 4: 295–306.
4. Tripodi A, Biasiolo A, Chantarangkul V, Pengo V. Lupus anticoagulant (LA) testing: performance of clinical laboratories assessed by a national survey using lyophilized affinity-purified immunoglobulin with LA activity. *Clin Chem* 2003; 49: 1608–14.
5. Jennings I, Kitchen S, Woods TA, Preston FE, Greaves M. Potentially clinically important inaccuracies in testing for the lupus anticoagulant: an analysis of results from three surveys of the UK National External Quality Assessment Scheme (NEQAS) for Blood Coagulation. *Thromb Haemost* 1997; 77: 934–7.
6. Favalaro EJ, Bonar R, Sioufi J, Wheeler M, Low J, Aboud M, Duncan E, Smith J, Exner T, Lloyd J, Marsden K. Multilaboratory testing of thrombophilia: current and past practice in Australasia as assessed through the Royal College of Pathologists of the Australasia Quality Assurance Programs for Hematology. *Semin Thromb Haemost* 2005; 31: 49–58.
7. Arnout J, Meijer P, Vermynen J. Lupus anticoagulant testing in Europe: an analysis of results from the first European Concerted Action on Thrombophilia (ECAT) survey using plasmas spiked with monoclonal antibodies against human beta2-glycoprotein I. *Thromb Haemost* 1999; 81: 929–34.
8. Pengo V, Biasiolo A, Gresele P, Marongiu F, Erba N, Veschi F, Ghirarduzzi A, de Candia E, Montaruli B, Testa S, Barcellona D, Tripodi A. on behalf of participating centers of Italian Federation of Thrombosis Centers (FCSA). Survey on lupus anticoagulant diagnosis by central evaluation of positive plasma samples. *J Thromb Haemost* 2007; 5: 925–30.

9. Clinical and Laboratory Standards Institute (CLSI). Collection, Transport, and Processing of Blood Specimens for Testing Plasma- Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document H21-A5 (ISBN 1-56238-657-)
10. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087–1898 USA, 2008. Favaloro EJ. Preanalytical variables in coagulation testing. *Blood Coagul Fibrinolysis* 2007; 18: 86–9.
11. Galli M, Finazzi G, Bevers EM, Barbui T. Kaolin clotting time and dilute Russell viper venom time distinguish between prothrombindependent and beta 2-glycoprotein I-dependent antiphospholipid antibodies. *Blood* 1995; 86: 617–23.
12. Urbanus RT, Derksen RH, de Groot PG. Current insight into diagnostics and pathophysiology of the antiphospholipid syndrome. *Blood Rev* 2008; 2: 93–105.
13. Rand JH, Wu XX, Quinn AS, Chen PP, Hathcock JJ, Taatjes DJ. Hydroxychloroquine directly reduces the binding of antiphospholipid antibody-beta2-glycoprotein I complexes to phospholipid bilayers. *Blood* 2008; 112: 1687–95.
14. Tripodi A, Chantarangkul V, Clerici M, Mannucci PM. Laboratory diagnosis of lupus anticoagulants for patients on oral anticoagulant treatment: performance of dilute Russell viper venom test and silica 88: 583–6.
15. Triplett DA, Barna LK, Unger GA. A hexagonal (II) phase phospholipids neutralization assay for lupus anticoagulant identification. *Thromb Haemost* 1993; 70: 787–93.
16. Jacobsen EM, Barna-Cler L, Taylor JM, Triplett DA, Wisloff F. The Lupus Ratio test: an interlaboratory study on the detection of Lupus anticoagulants by an APTT-based, integrated, and semiquantitative test. Fifth International Survey of Lupus Anticoagulants ISLA 5. *Thromb Haemost* 2000; 83: 704–8.
17. Rosner E, Pauzner R, Lusky A, Modan M, Many A. Detection and quantitative evaluation of lupus circulating anticoagulant activity. *Thromb Haemost* 1987; 57: 144–7.
18. Triplett DA, Stocker KF, Unger GA, Barna LK. The Textarin/Ecarin ratio: a confirmatory test for lupus anticoagulants. *Thromb Haemost* 1993; 70: 925–31.
19. Moore GW, Smith MP, Savidge GF. The Ecarin time is an improved confirmatory test for the Taipan snake venom time in warfarinized patients with lupus anticoagulants. *Blood Coagul Fibrinolysis* 2003; 14: 307–12.
20. Pengo V, Biasiolo A, Gresele P, Marongiu F, Erba N, Veschi F, Ghirarduzzi A, Barcellona D, Tripodi A. A comparison of lupus anticoagulant-positive patients with clinical picture of antiphospholipid syndrome and those without. *Arterioscler Thromb Vasc Biol* 2007; 27: 309–10. 1740 V. Pengo et al