



## REAGENT STORAGE

The reconstituted Lupus Anticoagulant Confirmation Reagent is stable for 30 days when stored at 2° - 8° C in its original tightly sealed container.

## TEST PROCEDURE

The platelet neutralization procedure (PNP) assay is based on the microsilica activator based activated partial thromboplastin time (APTT) technique. IT SHOULD BE NOTED THAT THE VARIABILITY IN APTT REAGENTS FROM DIFFERENT MANUFACTURERS CAN HAVE AN EFFECT ON THE PNP.<sup>16-22</sup> The APTT methodology currently in use in the laboratory should be followed.

### PART 1: ROUTINE BASELINE APTT

1. Pipette 1 part of normal control plasma into a test cuvette.
2. Add 1 part of APTT reagent to the control plasma. Mix well.
3. Incubate the plasma-APTT mixture for the activation time recommended by the reagent manufacturer.
4. Add 1 part of prewarmed 0.025 M calcium chloride, simultaneously starting the timer.
5. Record the clotting time.
6. Repeat steps 1-5 for a duplicate sample. Duplicate results should correlate within  $\pm 5\%$ .
7. Record results for final evaluation.
8. Repeat steps 1-7 for each test plasma.

### PART 2: PLATELET NEUTRALIZATION PROCEDURE

1. Pipette 1 part normal control plasma into a test cuvette.
2. As a dilution control, add 0.1mL TRIS Buffered Saline, 0.06 M pH 7.5 or 0.85% (w/v) saline.
3. Add 1 part APTT reagent, mix well.
4. Incubate for recommended activation time, as in Part 1, Step 3.
5. Add 1 part prewarmed calcium chloride, simultaneously starting the timer.
6. Repeat steps 1 - 5 for a duplicate sample. Duplicate results should correlate within  $\pm 5\%$ . Record saline/normal control plasma results for final evaluation.
7. Repeat steps 1 - 6 for each test plasma. Record saline control for each test plasma.
8. Pipette 1 part normal control plasma into test cuvette.
9. Add 1 part of reconstituted Lupus Anticoagulant Confirmation Reagent.
10. Add 1 part of APTT reagent. Mix well.
11. Incubate for recommended activation time.
12. Add 1 part prewarmed calcium chloride, simultaneously starting timer.
13. Repeat steps 8 - 12 for duplicate sample. Duplicate results should correlate within  $\pm 5\%$ . Record Lupus Anticoagulant Confirmation Reagent/normal control plasma results for final evaluation.
14. Repeat steps 8 -13 for each test plasma. Record each Lupus Anticoagulant Confirmation Reagent/ test plasma result for final evaluation.

## QUALITY CONTROL

Performance of platelet neutralization procedure on normal control plasma is necessary to provide a reference point for interpretation of patient test results. The control should be assayed exactly as the test plasma. APTT values should be ascertained on a 1:1 mixture of saline/normal control plasma and a 1:1 mixture of Lupus Anticoagulant Confirmation Reagent/ normal control plasma.

## RESULTS AND FINAL EVALUATION

The clotting time of the Lupus Anticoagulant Confirmation Reagent/test plasma mixture is compared to its saline/test plasma control. A correction of the prolonged APTT of 5 seconds or greater for the Lupus Anticoagulant Confirmation Reagent/test plasma mixture as compared to the saline/ test plasma control, is considered a positive test result.

## EXPECTED VALUES

Expected Ranges for coagulation testing should be established by each laboratory.

## CORRECTION DETERMINATION

Subtract the LACR APTT time (secs) from the saline APTT time (secs) to determine the correlation

Saline Test Plasma APTT - LACR Test Plasma = correction

An example of three correction study panels is illustrated:

Patient	Baseline APTT (sec)	APTT (sec) saline/test plasma	APTT (sec) LA Conf Rgt /test plasma	Difference between APTT/LA Conf Rgt & APTT/saline	Eval.
1	49	44	43	1	Neg.
2	50	53	52	1	Neg.
3	75	53	41	12	Pos.

Patients 1 and 2 are reported as negative for the lupus anticoagulant. Patient 1 results are indicative of a specific factor inhibitor rather than a lupus anticoagulant because the APTT diluted with Lupus Anticoagulant Confirmation Reagent and the APTT diluted with saline BOTH corrected from the baseline APTT to the same extent. Here, correction is attributed to dilution of an inhibitor. Patient 2 APTT test results showed no correction with the addition of the Lupus Anticoagulant Confirmation Reagent or saline. The lack of correction is interpreted as negative for the lupus anticoagulant and more indicative of the presence of a factor deficiency rather than an inhibitor. Patient 3 results are indicative of the presence of a lupus anticoagulant. Refer to the flow chart for instructions on the evaluation of a prolonged APTT.

It should be noted that expected values may vary from laboratory to laboratory due to the heterogeneity of laboratory results in patients with the lupus anticoagulant.<sup>16</sup> Additionally, differences in coagulation results have been described to occur due to sensitivity and responsiveness of activated partial thromboplastin time reagents.<sup>17,18,19,20</sup> Parameters for expected values should be established by each laboratory.

## LIMITATIONS

Lupus Anticoagulant Confirmation Reagent is a phosphatidyl enhanced platelet phospholipid used to determine if a prolonged APTT is caused by a lupus inhibitor or another coagulation factor deficiency.<sup>19</sup> Lupus anticoagulants demonstrate considerable heterogeneity and show variable differences in sensitivity and responsiveness to the APTT reagent. The sensitivity and responsiveness may vary from lot to lot of APTT reagent. The quality of a result can only be as good as the quality of the specimen. To avoid erroneous results, specimen should be free of platelets and hemolysis.

Clinical history is required for accurate test interpretation. Patients should be questioned about the recent ingestion of any medication because a number of prescription and non-prescription drugs may interfere with coagulation. Dose and time taken should be noted for patients on Coumadin®.

## PERFORMANCE CHARACTERISTICS

The Lupus Anticoagulant Confirmation Reagent has been tested on plasmas of patients with a lupus anticoagulant. Studies have shown that the shortening effect on the APTT may vary with the titer of the inhibitor in the test plasma.

Studies have shown that this product will perform as described when procedural and storage directions

are followed.

## Linearity

The complex chemistry and nature of measurement of Prothrombin Time (PT) and Activated Partial

Thromboplastin Time (APTT) testing is not a linear relationship. Clotting times do not change in a linear ratio to the clotting factor(s) that are being measured. The limits of the instrument and method dictate the linearity of this assay

## Accuracy, Precision and Reproducibility

The accuracy, precision and reproducibility of coagulation testing are technique and instrument dependent. The Laboratory should establish its own limits of acceptability based on written laboratory protocols and accepted laboratory standards.

## REFERENCES

1. Nygaard, KK. Hemorrhagic Diseases: Photo-Electric Study of Blood Coagulation. C.V. Mosby Company, St. Louis. 1941. p17-21.
2. Siegel JD, Rhinehart E, Jackson M, Chiarello L, and the Healthcare Infection Control Practices Advisory Committee, 2007 Guidelines for Isolation Precautions; Preventing Transmission of Infectious Agents in Healthcare Settings <http://cdc.gov/ncidod/dhqp/pdf/isolation2007/pdf>
3. Clinical and Laboratory Standards Institute (CLSI). Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline- Third Edition. CLSI document M29-A3 (ISBN 1-56238-567-4). Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
4. Clinical and Laboratory Standards Institute (CLSI). Collection, Transport and processing of blood Specimens for testing Plasma-Based Coagulation assays and Molecular Hemostasis Assays. Approved Guideline-Fifth Edition. CLSI document H21 -A5. (ISBN 1-56238-657-3). Clinical and Laboratory Standards Institute 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2008.
5. Triplett DA, Brandt JT, Kaczor D, Schaeffer JS: Laboratory diagnosis of lupus inhibitors: A comparison of the tissue thromboplastin inhibition procedure with a new platelet neutralization procedure. Am J Clin Path 1982; 79: 678-682.
6. Green D. et al: Report of the working party on acquired inhibitors of coagulation studies of the "lupus" anticoagulant: Thromb Haemost 1983; 49: 144-146.
7. Schleider M. et al: A clinical study of the lupus anticoagulant, Blood 1976; 48: 499-509.
8. Exner T, Richard KA, Kronenberg H: A sensitive test demonstrating lupus anticoagulant and its behavioral pattern. Br J Haematol 1978; 40: 143-151.
9. Gastineau DA, et al: Lupus anticoagulant: an analysis of the clinical laboratory features of 219 cases. Am J Hematol 1985; 19: 265-275.
10. Thaigaranjan P, Pengo V, Shapiro SS: The use of the Russell viper venom time for the diagnosis of lupus anticoagulant. Blood 1986; 68: 869-874.
11. Alving BM, et al: The dilute phospholipid APTT: a sensitive assay for verification of lupus anticoagulants. Thromb Haemost 1985; 54: 709-712.
12. Kelsey PR: The diagnosis of LA by the activated partial thromboplastin time: The central role of phosphatidyl serine. Throm Haemost 1984; 52: 172-175.
13. Espinoza LR: Significance of the lupus anticoagulant. Am J Hematol 1986; 22: 331-337.
14. Rosove MH, et al: lupus anticoagulants: Improved diagnosis with a kaolin clotting time using rabbit brain phospholipid in standard and high concentrations. Blood 1986; 68: 472-478.
15. Kelsey PR, Stevenson KJ: lupus anticoagulant and phosphatidyl serine: demonstration of the interrelationship and a possible simple diagnostic test. (Abstract) Br J Haematol 1984; 56: 671.
16. Triplett DA, Brandt JT, Maas RL: The laboratory heterogeneity of lupus anticoagulants. Arch Path Lab Med 1985; 109: 946-951.
17. Mannucci PM, Canciani MT, Mari D, et al: The varied sensitivity of partial thromboplastin and prothrombin reagents in the demonstration of the lupus-like anticoagulant. Scand J Haematol 1979; 22: 423-432.
18. Brandt JT, Triplett DA, Musgrave K, Orr C: The sensitivity of different coagulation reagents in the presence of lupus anticoagulants. Arch Patol Lab Med 1987; 111: 120-124.
19. Brandt JT, Triplett DA, Rock WA, Bovill EG, Arkin CF: Effect of lupus anticoagulants on the activated partial thromboplastin time: results of coll of Amer Path Survey Program. Arch Path and Lab Med 1991; 115: 109-114.
20. Malia RG, Cooper SM, Hill FE, Preston M: Identification of lupus anticoagulant in patients receiving oral anticoagulants. Abstract XIII Congress of Int Soc on Thrombosis and Hemostasis. Amsterdam June 1991.
21. Triplett DA, Hannah D, Barna L: Comparison of platelet extracts for in-vitro neutralization of lupus anticoagulant. Abstract XIII Congress of Int Soc on Thrombosis and Hemostasis. Amsterdam June 1991.
22. Brian JT, Shaus K, Cooper K, Shaus M, O'Keefe B: Evaluation of platelet neutralization procedure using different APTT reagents. Abstract XIII Congress of Int Soc on Thrombosis and Hemostasis. Amsterdam June 1991.
23. Harmening D. Clinical Hematology and Fundamentals of Hemostasis, 3rd edition, 543-546, FA Davis Company, Phila, 1997.
24. Diagnosis of Hemophilia and Other Bleeding Disorder: A LABORATORY MANUAL, Second Edition. World Federation of Hemophilia, Montreal, Quebec, CANADA (2010).
25. Pengo, V, Tripodi, A, Reber, G, Rand, JH, Ortel, TL, Galli, M, and DeGroot, PG. Update of the Guidelines for Lupus Anticoagulant Detection. Official Communication of the SSC. J Thromb Haemost, 7: 1737 - 1740. 2009.
26. Utz, VM and Tang, J. Ocular Manifestations of the Antiphospholipid Syndrome. Br J Ophthalmol. Doc 10.1136/bjo.182857, 2010
27. Hill, GS and Nochy, D. Antiphospholipid Syndrome in Systemic Lupus Erythematosus. J AM Soc Neph.18: 2461 - 2464. 2007

For a complete list of available products please go to [www.biodatacorp.com](http://www.biodatacorp.com)

THIS PRODUCT IS WARRANTED TO PERFORM AS DESCRIBED IN THE LABELING AND IN THE LITERATURE OF BIO/DATA CORPORATION AND BIO/DATA CORPORATION DISCLAIMS ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY OTHER PURPOSE AND IN NO EVENT SHALL BIO/DATA CORPORATION BE LIABLE FOR ANY CONSEQUENTIAL DAMAGES ARISING OUT OF THE AFORESAID EXPRESSED WARRANTY.



155 Gibraltar Road, Horsham, PA 19044 U.S.A.  
(800) 257-3282 U.S.A. (215) 441-4000 Worldwide  
(215) 443-8820 Fax Worldwide  
Email: [customer.service@biodatacorp.com](mailto:customer.service@biodatacorp.com)  
Internet: [www.biodatacorp.com](http://www.biodatacorp.com)



EMERGO EUROPE, Prinsessegracht 20, 2514 AP, The Hague, The Netherlands

