

# The PAP-8E Quick Guide

## GENERAL TEST PROCEDURE FOR THE PAP-8E AGGREGOMETER

1. Prepare a blank by pipetting 250 $\mu$ L of platelet poor plasma (PPP) into a test tube - do not add a stir bar
  - a) Place blank in a test tube rack
2. Place an empty test tube into an incubation well - add a stir bar and wait at least one minute
3. Pipette 225 $\mu$ L of platelet rich plasma (PRP) into the pre-warmed test tube
4. After adding the PRP to the pre-warmed test tube
  - a) Select **Timer** and the count down will start
  - b) Allow the sample to incubate at 37° C for the preset time
5. Set the 100% baseline by placing the blank into the test well
  - a) Select **Blank**
  - b) When the Blank is set, the Blank button will change to Start
6. Place the tube of PRP into the test well and close the pipette guide
  - a) Select **Start**
7. Add 25 $\mu$ L of reagent directly into the PRP
  - a) Select **Inject**
8. The test will run for a pre-set time

For Additional Procedures or Directions,  
Click **Help** For The Onboard Manual



155 Gibraltar Road  
Horsham, PA 19044 USA

Worldwide: +1 215-441-4000  
USA: 1-800-257-3282  
Fax Worldwide: +1 215-443-8820  
customer.service@biodatacorp.com

REF  
106075  
106077



An ISO 13485 Registered Company  
www.biodatacorp.com

Proudly Manufactured in the USA



mdi Europa GmbH  
Langenhagener Str. 71  
D-30855 Langenhagen  
GERMANY

Alpha Laboratories  
40 Parham Drive Eastleigh  
SO50 4NU Hampshire  
UNITED KINGDOM

