

PRODUCT DESCRIPTION

LTA Check™ ✓ is a system of reagents designed to verify the proper operation and performance of Light Transmission Aggregometers.

LTA Check™ ✓ Performance Monitoring Test Kit has been optimized for use with Light Transmission Aggregometers.

INTENDED PURPOSE

LTA Check™ ✓ is a Performance Monitoring Test Kit used to verify the operation of Light Transmission Aggregometers (LTA) and provide instrument system quality control in the laboratory.

DETECTION / MEASUREMENT

LTA Check™ ✓ reagents are used, in conjunction with the provided diluent, to measure changes in light transmission within a controlled test system.

PRODUCT FUNCTION

LTA Check™ ✓ Performance Monitoring Test Kit is designed to verify the performance and proper operation of Light Transmission Aggregometers. The system generates defined aggregation and slope responses that enable evaluation of instrument functionality, precision, and reproducibility under controlled conditions.

SPECIFIC INFORMATION PROVIDED

LTA Check™ ✓ Performance Monitoring Test Kit is not intended for the detection of a specific disorder, condition, or risk factor.

This test kit provides information related to the performance of Light Transmission Aggregometers, including the consistency, precision, and reproducibility of aggregation parameters generated under defined test conditions. The results reflect instrument system functionality and are used to verify proper operation of the aggregation system.

AUTOMATION

LTA Check™ ✓ Performance Monitoring Test Kit is intended for use with semi-automated and automated Light Transmission Aggregometers (LTA).

QUALITY / QUANTITY

There are no primary standards for LTA Check™ ✓ Performance Monitoring Test Kit. The responses generated by this system are concentration dependent.

LTA Check™ ✓ Performance Monitoring Test Kit comes packaged as 1 x 0.5 mL vial of Reagent 1, 1 x 4.0 mL vial of Reagent 2, 1 x 1.0 mL vial of Reagent 3, 1 x 10.0 mL vial of LTA Diluent, and 1 x 1.0 mL vial of Ultra Pure Water.

SPECIMEN TYPE

No patient specimen is required for use with the LTA Check™ ✓ Performance Monitoring Test Kit. The test system utilizes the reagents and diluent provided within the kit to generate controlled reactions for performance verification of Light Transmission Aggregometers.

TESTING POPULATION

LTA Check™ ✓ Performance Monitoring Test Kit is not intended for use with a patient population. The test kit is intended for use by trained laboratory personnel to verify the performance and proper operation of Light Transmission Aggregometers within professional laboratory settings.

IN VITRO DIAGNOSTIC

LTA Check™ ✓ Performance Monitoring Test Kit is an in vitro diagnostic test kit intended for Professional Laboratory Use Only. It is not intended for injection or ingestion.

INTENDED USER

LTA Check™ ✓ Performance Monitoring Test Kit is intended for Professional Laboratory Use by qualified personnel.

TEST PRINCIPLE

LTA Check™ ✓ Performance Monitoring Test Kit evaluates the performance of Light Transmission Aggregometers by generating a controlled agglutination reaction within the test system. When the reagents are combined under defined conditions and stirred at 37°C, agglutination occurs, resulting in measurable changes in light transmission. The Light Transmission Aggregometer detects these changes and expresses the reaction as functional parameters, including Final Aggregation (FA) and Slope (PS), which represent the extent and rate of the reaction. These parameters are used to verify the proper operation and performance of the instrument system.

CALIBRATORS AND CONTROLS

There are no external calibrators required for the LTA Check™ ✓ Performance Monitoring Test Kit. The reagents provided within the kit function as a performance monitoring system and are used to generate defined aggregation parameters for evaluation of instrument performance. Each laboratory should establish its own expected ranges and acceptance criteria for aggregation parameters when using this test kit.

REAGENT LIMITATIONS

LTA Check™ ✓ Performance Monitoring Test Kit will perform as specified when the Instructions for Use are followed. The reagents must be used prior to the expiration date printed on each vial.

REAGENTS PROVIDED

REF	107117:	1 vial of Reagent 1 (0.5 mL)
		1 vial of Reagent 2 (4.0 mL)
		1 vial of Reagent 3 (1.0 mL)
		1 vial of LTA Diluent (10.0 mL)
		1 vial of Ultra Pure Water (1.0 mL)

REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

No additional reagents are required for use with the LTA Check™ ✓ Performance Monitoring Test Kit. All necessary reagents and diluents are provided.

MATERIALS AND ACCESSORIES

- Platelet Aggregometer (*Follow the Manufacturer's Instructions for Use*)
- Electronic Pipette
- Pipette Tips ②
- Aggregometer Test Tubes (*Siliconized*) ②
- Aggregometer Stir Bars (*Plastic Coated*) ②
- Plastic Sample Tubes and Caps (*for Dilutions*) ②



NOTE: DISPOSABLE ITEMS SUCH AS TEST TUBES, STIR BARS, SAMPLE TUBES, AND CAPS ARE FOR ONE TIME USE ONLY







STORAGE AND STABILITY

- LTA Check™ ✓ Performance Monitoring Test Kit does not require temperature protection during shipment.
- Upon receipt, store Reagent 1, Reagent 2, Reagent 3, LTA Diluent, and Ultra Pure Water at 2 – 8° C in their original packaging.
- Reconstituted Reagent 1 is stable for 7 days when stored in its tightly capped, original container at 2 – 8° C.
- Reconstituted Reagent 2 is stable for 30 days when stored in its tightly capped, original container at 2 – 8° C.
- Reconstituted Reagent 3 is stable for 8 hours when stored in its tightly capped, original container at 2 – 8° C.

STERILITY

LTA Check™ ✓ Performance Monitoring Test Kit is not a sterile product. Be careful not to contaminate the product when pipetting the reconstituted or aliquoted reagents.

WARNINGS AND PRECAUTIONS

-  Wear PPE in accordance with laboratory policies and practices when handling LTA Check™ ✓ Performance Monitoring Test Kit.
-  Follow standard precautions when preparing test specimens and samples.
-  Handle LTA Check™ ✓ Performance Monitoring Test Kit with care to avoid contamination during use.
-  Avoid reagent evaporation by limiting air – liquid exchange surfaces.
-  To preserve reagent stability, store remaining reagent in its tightly capped, original container.
-  Dispose of post-test materials in accordance with applicable regulations and laboratory policies.



NOTE TO USER: ANY SERIOUS INCIDENT THAT OCCURS IN RELATION TO THIS PRODUCT SHALL BE REPORTED TO THE MANUFACTURER AND THE COMPETENT AUTHORITY OF THE MEMBER STATE IN WHICH THE USER AND / OR PATIENT ARE ESTABLISHED.

INFECTIOUS MATERIAL STATUS

LTA Check™ ✓ Reagent 1, LTA Diluent, and Ultra Pure Water do not contain any infectious materials. Reagent 2 and Reagent 3 contain plasma or platelets that have been screened at the source and confirmed negative for HIV-1 antigen (HIV-1Ag), anti-HIV-1/2 antibodies, Hepatitis B surface antigen (HBsAg), Hepatitis C antibodies, Human T-Lymphotropic Virus Type I and II (anti-HTLV I/II) antibodies, and Syphilis by serological testing. Nevertheless, due to their human origin, all plasma and platelets should be handled as potentially hazardous materials. After testing, test specimens and samples must be disposed of in compliance with applicable regulations and laboratory policies.

SPECIAL FACILITIES

LTA Check™ ✓ Performance Monitoring Test Kit does not require the use of special facilities within a laboratory environment.

PREPARATION FOR USE



NOTE: LTA CHECK™ ✓ PERFORMANCE MONITORING TEST KIT MUST BE AT ROOM TEMPERATURE (15 – 28° C) PRIOR TO RECONSTITUTION. STORED REAGENTS MUST BE BROUGHT TO ROOM TEMPERATURE PRIOR TO USE.

RECONSTITUTION

- Reconstitute Reagent 1 with 0.5 mL of Ultra Pure Water.
- Reconstitute Reagent 2 with 4.0 mL of LTA Diluent.
- Reconstitute Reagent 3 with 1.0 mL of LTA Diluent.
- Cap each vial.
- Wait 5 minutes after reconstituting each vial.
- Invert each vial gently to mix.
- Allow each vial to stand at room temperature for 25 minutes.
- Invert each vial gently to ensure complete mixing and rehydration.
- Reconstituted vials should be kept capped prior to use.



NOTE: TESTING MUST BE COMPLETED WITHIN 45 MINUTES OF REAGENT PREPARATION.

PATIENT PREPARATION

LTA Check™ ✓ Performance Monitoring Test Kit does not require patient preparation, as no patient specimens are used.

SPECIMEN COLLECTION

LTA Check™ ✓ Performance Monitoring Test Kit does not require specimen collection, as no patient specimens are used.



PRACTICE STANDARD PRECAUTIONS THROUGHOUT THE SPECIMEN COLLECTION, SAMPLE PREPARATION, AND ANALYTICAL PROCESSES. DISPOSE OF SHARPS AND BIOHAZARDOUS WASTE IN ACCORDANCE WITH APPLICABLE REGULATIONS AND LABORATORY POLICIES.

SAMPLE PREPARATION

LTA Check™ ✓ Performance Monitoring Test Kit does not require preparation of patient specimens. Reagent preparation is described in the RECONSTITUTION section.

ASSAY PROCEDURE



NOTE: THIS IS A GENERAL PROCEDURE. FOLLOW THE INSTRUCTIONS FOR USE PROVIDED BY THE MANUFACTURER OF THE AGGREGOMETER IN USE.



NOTE: TESTING MUST BE COMPLETED WITHIN 45 MINUTES OF REAGENT PREPARATION.

FULL VOLUME PROCEDURE 500 µL

Prepare a Blank

- Label a test tube with the letter "B" to identify the Blank
- Pipette 250 µL of Reagent 2 and 250 µL of LTA Diluent into the test tube (DO NOT ADD A STIR BAR)
- Invert gently to mix
- Let stand for 10 minutes
- Invert gently again to mix prior to use

Prepare Samples

- Label one to eight new test tubes with test well #
- Place the labeled test tubes into the correct well # 1 - 8 of the stirred sample incubation wells
- Add a stir bar to each test tube
- Incubate test tubes for one minute
- Pipette 50 µL of reconstituted Reagent 1 into each test tube in the stirred sample incubation wells (MAKE SURE THERE ARE NO BUBBLES)
- Pipette 400 µL of reconstituted Reagent 2 into each test tube in the stirred sample incubation wells (MAKE SURE THERE ARE NO BUBBLES)
- Select the on-screen timer for each stirred sample incubation well in use and

the warming countdown will start

- The samples will incubate at 37° C for the pre-set time (2 Minutes)

Set the 100% baseline (Blank)

- Place the previously prepared Blank test tube into test well # 1
- Select BLANK to activate the test well
- The BLANK button will change to START
- Repeat the steps above for each test well being used for testing

Begin Testing

- Once the countdown timer reaches 0:00, press the timer button to stop each stirred sample incubation well
- Transfer the test tube in the stirred sample incubation well # 1 to test well # 1
- Repeat the step above for each test well, making sure all test tubes remain with their corresponding well #'s during transfer
- Close the pipette guides
- Select START for test well # 1
- Pipette 50 µL of Reagent 3 directly into the test tube in test well # 1 (DO NOT ALLOW REAGENT TO RUN DOWN THE WALL OF THE TEST TUBE OR PERMIT PIPETTE TIP TO BREAK THE SURFACE OF THE SAMPLE)
- Select INJECT for test well # 1
- Repeat the steps above for each test well being used for testing
- The test will now run for the pre-set time (5 Minutes) (OTHER MANUFACTURER'S TEST PROCEDURES MAY SPECIFY DIFFERENT TIMES OR VOLUMES)



NOTE: EACH LABORATORY SHOULD ESTABLISH AND VALIDATE ITS OWN TEST PROTOCOL AND VERIFY THE RESULTING PERFORMANCE OF ITS TEST SYSTEM (REAGENTS, INSTRUMENT, AND TEST PROTOCOL).

HALF VOLUME PROCEDURE 250 µL

Prepare a Blank

- Label a test tube with the letter "B" to identify the Blank
- Pipette 125 µL of Reagent 2 and 125 µL of LTA Diluent into the test tube (DO NOT ADD A STIR BAR)
- Invert gently to mix
- Let stand for 10 minutes
- Invert gently again to mix prior to use

Prepare Samples

- Label one to eight new test tubes with test well #
- Place the labeled test tubes into the correct well # 1 - 8 of the stirred sample incubation wells
- Add a stir bar to each test tube
- Incubate test tubes for one minute
- Pipette 25 µL of reconstituted Reagent 1 into each test tube in the stirred sample incubation wells (MAKE SURE THERE ARE NO BUBBLES)
- Pipette 200 µL of reconstituted Reagent 2 into each test tube in the stirred sample incubation wells (MAKE SURE THERE ARE NO BUBBLES)
- Select the on-screen timer for each stirred sample incubation well in use and the warming countdown will start
- The samples will incubate at 37° C for the pre-set time (2 Minutes)

Set the 100% baseline (Blank)

- Place the previously prepared Blank test tube into test well # 1
- Select BLANK to activate the test well
- The BLANK button will change to START
- Repeat the steps above for each test well being used for testing

Begin Testing

- Once the countdown timer reaches 0:00, press the timer button to stop each stirred sample incubation well
- Transfer the test tube in the stirred sample incubation well # 1 to test well # 1
- Repeat the step above for each test well, making sure all test tubes remain with their corresponding well #'s during transfer
- Close the pipette guides
- Select START for test well # 1
- Pipette 25 µL of Reagent 3 directly into the test tube in test well # 1 (DO NOT ALLOW REAGENT TO RUN DOWN THE WALL OF THE TEST TUBE OR PERMIT PIPETTE TIP TO BREAK THE SURFACE OF THE SAMPLE)
- Select INJECT for test well # 1
- Repeat the steps above for each test well being used for testing
- The test will now run for the pre-set time (5 Minutes) (OTHER MANUFACTURER'S TEST PROCEDURES MAY SPECIFY DIFFERENT TIMES OR VOLUMES)



NOTE: EACH LABORATORY SHOULD ESTABLISH AND VALIDATE ITS OWN TEST PROTOCOL AND VERIFY THE RESULTING PERFORMANCE OF ITS TEST SYSTEM (REAGENTS, INSTRUMENT, AND TEST PROTOCOL).

QUALITY CONTROL

Quality control is achieved by generating defined aggregation parameters, including Final Aggregation (FA) and Primary Slope (PS), and evaluating the precision and reproducibility of these results under consistent test conditions.

Testing should be performed periodically, with new reagent lots, and following instrument maintenance to verify consistent system performance. Each laboratory should

establish its own expected ranges and acceptance criteria for aggregation parameters and verify ongoing system performance.

RESULTS

Results are expressed as functional parameters generated by the Light Transmission Aggregometer, including Final Aggregation (FA) and Primary Slope (PS), which represent the extent and rate of the reaction.

Results should be evaluated based on the precision and reproducibility of these parameters under defined test conditions. Consistent results indicate proper instrument performance, while variability may indicate issues with instrument function, reagent handling, or test conditions.

Each laboratory should establish its own expected ranges and acceptance criteria for aggregation parameters and use these to assess system performance.

LIMITATIONS

The performance of the LTA Check™ ✓ Performance Monitoring Test Kit is dependent on proper instrument operation, reagent handling, and adherence to specified test conditions.

Variability may occur due to differences in instrument configuration, environmental conditions, and operator technique.

This test kit is intended for performance verification of Light Transmission Aggregometers and does not provide information related to patient conditions.

EXPECTED VALUES

Expected results should be evaluated based on the precision and reproducibility of aggregation parameters. Typical performance characteristics are shown below.

Expected CV (vs. average of established range):

- Expected Primary Slope: 25–50 (± 10%)
- Expected Final Aggregation: 60–110% (± 10%)

These values were established using Bio/Data Corporation's Platelet Aggregation Profiler, Model PAP-8E. Other instruments or models may yield different results. Each laboratory should establish its own expected ranges and acceptance criteria.

ANALYTICAL PERFORMANCE

LTA Check™ ✓ Performance Monitoring Test Kit is designed to evaluate the functionality of Light Transmission Aggregometers under controlled conditions. The test system generates aggregation responses that are measured as functional parameters, including Final Aggregation (FA) and Primary Slope (PS).

The performance of the system is dependent on proper instrument operation, reagent handling, and adherence to specified test conditions. Variability may be introduced by factors such as instrument configuration, environmental conditions, and operator technique.

Precision and reproducibility are assessed through the consistency of aggregation parameters generated during repeated testing. The following performance characteristics are typical for the system:

Test to Test Reproducibility:	less than ± 7.5%
Instrument to Instrument Reproducibility:	less than ± 15.0%
Reagent Lot to Lot Variability:	less than ± 10.5%
Laboratory to Laboratory (System to System)	less than ± 12.5%

OPTIONAL SOFTWARE

LTA Check™ ✓ Performance Monitoring Software (C/N 107149) is available for use with this test kit. The software accepts data generated by the Light Transmission Aggregometer and provides automated analysis, including the generation of printable Levey-Jennings charts for evaluation and review of system performance.

Levey-Jennings charts provide a graphical representation of results over time, allowing for visual assessment of whether performance remains within laboratory-established limits. The software supports the application of Westgard rules to enhance error detection and reduce false rejection of valid results.

REFERENCES

- Allain JP, Cooper HA, Wagner RH, Brinkhous KM. Platelets fixed with paraformaldehyde: a new reagent for assay of von Willebrand factor and platelet aggregating factor. *J Lab Clin Med.* 1975 Feb;85(2):318-28.
- Angiolillo DJ, Ueno M, Goto S. Basic principles of platelet biology and clinical implications. *Circ J.* 2010 Apr;74(4):597-607.
- Born GV, Cross MJ. The Aggregation of Blood Platelets. *J Physiol.* 1963 Aug; 168(1):178-95.
- Brinkhous KM, Read MS. Preservation of platelet receptors for platelet aggregating factor/von Willebrand factor by air drying, freezing, or lyophilization: new stable platelet preparations for von Willebrand factor assays. *Thromb Res.* 1978 Oct;13(4):591-7.
- Bye A, Lewis Y, O'Grady J. Effect of a single oral dose of aspirin on the platelet

aggregation response to arachidonic acid. *Br J Clin Pharmacol.* 1979 Mar; 7(3):283-6.

- Cattaneo M, Cerletti C, Harrison P, Hayward CP, Kenny D, Nugent D, Nurden P, Rao AK, Schmaier AH, Watson SP, Lussana F, Pugliano MT, Michelson AD. Recommendations for the Standardization of Light Transmission Aggregometry: A Consensus of the Working Party from the Platelet Physiology Subcommittee of SSC/ISTH. *J Thromb Haemost.* 2013 Apr 10.
- CLSI. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. CLSI document H18-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
- CLSI. Protection of Laboratory Workers from Occupationally Acquired Infections, Approved Guideline - Fourth Edition. CLSI document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
- CLSI. Platelet Function Testing by Aggregometry, Approved Guideline - Fourth Edition. CLSI document H58-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- CLSI. Collection, Transport and Processing for Plasma Based Coagulation Assays and Molecular Hemostasis Assays, Approved Guideline - Fifth Edition. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- CLSI. Clinical Laboratory Safety, Approved Guideline - Third Edition. CLSI document GP17-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- Day HJ, Holmsen H. Laboratory tests of platelet function. *Ann Clin Lab Sci (1971).* 1972 Jan-Feb; 2(1):63-74.
- Day HJ, Rao AK. Evaluation of platelet function. *Semin Hematol.* 1986 Apr;23(2):89-101.
- Eichelberger, JW. Kinetic (Slope) Measurement of Platelet Aggregation. Bio/Data Corporation, Horsham, PA; 1984.
- Favalaro EJ, Gosselin RC, Pasalic L, Lippi G. Post-analytical issues in hemostasis and thrombosis testing: An update. In EJF, RCG, editors, Hemostasis and Thrombosis: Methods and Protocols. 2nd ed. New York: Humana Press. 2023. p. 787-811. (Methods in Molecular Biology).
- Federici AB, Lee CA, Berntorp EE, Lillicrap D, Montgomery RR. Von Willebrand Disease: Basic and Clinical Aspects. 2011.
- Garner JS. Guideline for isolation precautions in hospitals. The Hospital Infection Control Practices Advisory Committee. *Infect Control Hosp Epidemiol.* 1996 Jan;17(1):53-80.
- Howard MA, Firkin BG. Ristocetin—a new tool in the investigation of platelet aggregation. *Thromb Diath Haemorrh.* 1971 Oct 31; 26(2): 362-9.
- Israels SJ, El-Ekiaby M, Quiroga T, Mezzano D. Inherited disorders of platelet function and challenges to diagnosis of mucocutaneous bleeding. *Haemophilia.* 2010 Jul;16 Suppl 5:152-9.
- Kambayashi J, Shinoki N, Nakamura T, Ariyoshi H, Kawasaki T, Sakon M, Monden M. Prevalence of impaired responsiveness to epinephrine in platelets among Japanese. *Thromb Res.* 1996 Jan 1;81(1):85-90.
- Levine PH. The effect of thrombocytopenia on the determination of platelet aggregation. *Am J Clin Pathol.* 1976 Jan;65(1):79-82
- Linnemann B, Schwonberg J, Mani H, Prochnow S, Lindhoff-Last E. Standardization of light transmittance aggregometry for monitoring antiplatelet therapy: an adjustment for platelet count is not necessary. *J Thromb Haemost.* 2008 Apr;6(4):677-83.
- Marcus AJ, Coleman RW, Hirsh J, Ivarer VJ, Salzman EW. Hemostasis and thrombosis: Basic Principles and Clinical Practice. Vol. 472. Philadelphia: JB Lippincott Company; 1982.
- McRedmond J. Platelet Aggregation Quality Control in a Multicenter, Phase 2 Clinical Trial. *Drug Information Journal.* 2011;45(6):731-740.
- Michelson, AD. Platelets. Third Edition. Amsterdam: Academic Press; 2013.
- Mills DC, Robb IA, Roberts GC. The release of nucleotides, 5-hydroxytryptamine and enzymes from human blood platelets during aggregation. *J Physiol.* 1968 Apr;195(3):715-29.
- Moncada S, Vane JR. Arachidonic acid metabolites and the interactions between platelets and blood-vessel walls. *N Engl J Med.* 1979 May 17;300(20):1142-7.
- NCCLS. Assays of von Willebrand Factor Antigen and Ristocetin Cofactor Activity; Approved Guideline. NCCLS document H51-A. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2002.
- O'Donnell CJ, Larson MG, Feng D, Sutherland PA, Lindpaintner K, Myers RH, D'Agostino RA, Levy D, Toeller GH; Framingham Heart Study. Genetic and environmental contributions to platelet aggregation: the Framingham heart study. *Circulation.* 2001 Jun 26;103(25):3051-6.
- Owen CA Jr, Bowie EJW, Thompson JH Jr. The Diagnosis of Bleeding Disorders. 2nd ed. Little, Brown, and Company; 1975.
- Palma-Barqueros V, Revilla N, Sánchez A, Zamora Cánovas A, Rodríguez-Alén A, Marín-Quílez A, González-Porras JR, Vicente V, Lozano ML, Bastida JM, Rivera J. Inherited Platelet Disorders: An Updated Overview. *Int J Mol Sci.* 2021 Apr 26;22(9):4521.
- Siegel JD, Rhinehart E, Jackson M, Chiarello L; Health Care Infection Control Practices Advisory Committee. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Health Care Settings. *Am J Infect Control.* 2007 Dec;35(10 Suppl 2):S65-164.
- The Hospital Infection Control Practices Advisory Committee, Centers for disease Control and Prevention, Public Health Service, U.S. Department of Health and Human Services. Guideline for isolation precautions in hospitals Part II. Recommendations for isolation precautions in hospitals. *American Journal of Infection Control.* 1996; Vol 24, Issue 1: 32-52.
- Triplett DA, et al. Platelet function: laboratory evaluation and clinical application. Chicago, IL: American Society for Clinical Pathology 1978.
- Weiss HJ. Aspirin and Platelets in Drugs and Hematologic Reactions. New York, NY: Dimitov and Nodine, eds. Grune and Stratton. 1974.
- Westgard JO, Carey RN, Wold S. Criteria for judging precision and accuracy in method development and evaluation. *Clin Chem.* 1974;20:825-833.
- Westgard JO, deVos DJ, Hunt MR, Quam EF, Carey RN, Garber CC. Method evaluation. American Society for Medical Technology, Bellaire, TX; 1978.

- Westgard JO, Groth T, Aronsson T, Falk H, deVerdier C-H. Performance characteristics of rules for internal quality control: probabilities for false rejection and error detection. Clin Chem. 1977;23:1857-1867.
- Westgard JO, Hunt MR. Use and interpretation of common statistical tests in method comparison studies. Clin Chem. 1973;19:49-57.
- White, M.M., and Jennings, L.K. Platelet Protocols: Research and Clinical Laboratory Procedures, Academic Press, Inc.; 1999.
- Williams WJ, Beutler E, Erslev AJ, Rundles RW. Hematology. New York, NY: McGraw-Hill. 1977.

REVISION HISTORY

Document No: 107135 Revision: AA, May 2026

- Modified Testing Instructions
- Implemented IVDR Regulatory Requirements
- Reformatted and Reconfigured to Enhance Operator Use

For a complete product catalog, please visit our website at www.biodatacorp.com or contact our Customer Service Department.

THE BIO/DATA CORPORATION PRODUCT LINE INCLUDES GENERAL PURPOSE, PROFESSIONAL LABORATORY USE REAGENTS INTENDED TO INDUCE AND REPORT PLATELET FUNCTION ACTIVITY AND RESPONSES. THIS PRODUCT IS WARRANTED TO PERFORM AS DESCRIBED IN ITS LABELING INCLUDING THE INSTRUCTIONS FOR USE. BIO/DATA CORPORATION MAKES NO CLAIM OR WARRANTY, EXPRESSED OR IMPLIED, OF THE CAPABILITY, FITNESS, OR MERCHANTABILITY FOR ANY OTHER PURPOSE. IN NO EVENT SHALL BIO/DATA CORPORATION BE LIABLE FOR ANY CONSEQUENTIAL DAMAGES ARISING OUT OF AFORESAID EXPRESSED WARRANTY.

SYMBOLS



Bio-Hazardous



Catalog Number



Caution



CE Marked & Registered Product



Consult Instructions For Use



European Union Representative



In Vitro Diagnostic Device



Manufacturer



Must Read



Non-Sterile



Single Use Only



Temperature Limitations



United Kingdom Marked & Registered Product



United Kingdom Representative

 155 Gibraltar Road
Horsham, PA 19044 USA

Worldwide: +1 215-441-4000
USA: 1-800-257-3282
FAX Worldwide: +1 215-443-8820
customer.service@biodatacorp.com

©BIO/DATA CORPORATION 2026



AN ISO 13485 REGISTERED COMPANY

www.biodatacorp.com

PROUDLY MANUFACTURED IN THE USA



mdi Europa GmbH
Langenhagener Str. 71
D-30855 Langenhagen GERMANY



Alpha Laboratories
40 Parham Drive Eastleigh
SO50 4NU Hampshire UNITED KINGDOM



LTA CHECK INSTRUCTIONS FOR USE # 107135 REV AA