

PRODUCT DESCRIPTION

vW Select™ is a Ristocetin Cofactor Activity Assay Test System containing AggRecetin Reagent (Ristocetin A Sulfate), AggRecetin Diluent, Lyophilized Platelets, vW Abnormal Control Plasma, vW Normal Control Plasma, vW Normal Reference Plasma, and TRIS Buffered Saline (TBS). The materials provided in vW Select™ are carefully chosen to enhance performance and ensure greater repeatability compared to standard component combinations used in traditional Ristocetin Cofactor Assay Test Kits.

AggRecetin Reagent is a lyophilized preparation of Ristocetin A Sulfate, a glycopeptide of unknown chemical structure which has been isolated from *Nocardia lurida*. Ristocetin contains in excess of 90% Ristocetin A.

AggRecetin Diluent is a specially formulated solution used to reconstitute lyophilized AggRecetin Reagent for Platelet Aggregation and Ristocetin Cofactor Activity Testing.

Lyophilized Platelets are a standardized and fixed suspension of human platelets, specifically designed for platelet aggregation testing. Used routinely in Ristocetin Cofactor Assay Activity Tests, these platelets retain the critical Gp1b receptor function, ensuring accurate and reliable results. The suspension of formalin-fixed human platelets simplifies the detection of von Willebrand Factor (vWF) inhibition in mixed patient and pooled plasma samples. The lyophilization process extends the shelf life.

vW Abnormal Control Plasma is lyophilized citrated human plasma which has been selectively depleted of von Willebrand Factor.

vW Normal Control Plasma is lyophilized citrated human plasma, confirmed to have a von Willebrand Factor activity level between 80% and 120%.

vW Normal Reference Plasma is lyophilized citrated human plasma that has tested negative for von Willebrand Factor (vWF) and has been standardized using World Health Organization (WHO) reference material. It is designed to provide a consistent von Willebrand Factor activity level between 90% and 110%. The vW Normal Reference Plasma Assay Sheet is provided, detailing the established Reference Range and Assay Value. Reference curve linearity and acceptability should be evaluated based on the supplied assay value.

TRIS Buffered Saline (TBS) is a prepared buffered saline 0.06M containing 0.01% Sodium Azide as a preservative.

vW Select™ has been optimized for use with Light Transmission Aggregometers. It may also be used with other turbidometric or impedance analyzers, and flow cytometers.

INTENDED PURPOSE

vW Select™ is an all-in-one Ristocetin Cofactor Activity Assay Test System designed for precise von Willebrand Factor assessment. This Test System includes AggRecetin Reagent, AggRecetin Diluent, Lyophilized Platelets, vW Abnormal Control Plasma, vW Normal Control Plasma, vW Normal Reference Plasma, TRIS-Buffered Saline (TBS), and lot-specific technical information. It is optimized to induce the agglutination of lyophilized platelets in the presence of the patient's Platelet Poor Plasma and Ristocetin.

DETECTION / MEASUREMENT

vW Select™ is used, in conjunction with other reagents, diluents and control samples, to measure changes of the light transmission in a Platelet Rich Plasma (PRP) test sample.

PRODUCT FUNCTION

vW Select™ provides essential insights into the functional activity of von Willebrand Factor (vWF). This Test System is designed to aid in the evaluation of von Willebrand Syndrome (vWS), supporting the diagnosis and classification of inherited or acquired vWF-related disorders.

SPECIFIC INFORMATION PROVIDED

vW Select™ is not intended for the detection of a specific disorder, condition, or risk factor.

AggRecetin Reagent plays a pivotal role in assessing the functional / quality activity of von Willebrand Factor (vWF) through the vW Cofactor Assay, a critical test for evaluating vWF activity and function. This assay relies on Ristocetin to induce the interaction between vWF and platelets, under shear stress conditions. By binding to vWF, Ristocetin triggers a conformational change in the protein, exposing binding sites for platelet glycoprotein Ib (GPIb), facilitating platelet adhesion and aggregation. The Ristocetin Cofactor activity of vWF, observed in the presence of Ristocetin Reagent, evaluates its ability to promote the agglutination of formalin-fixed or fresh platelets in a patient's plasma. A decrease in Ristocetin Cofactor activity indicates an impairment or dysfunction of vWF, characteristic of von Willebrand Disease (vWD) or related disorders, by assessing the extent of platelet agglutination in response to Ristocetin.

AggRecetin Diluent is optimized and quality-controlled to ensure compatibility with AggRecetin Reagent included in vW Select™. It is used exclusively for the reconstitution of AggRecetin Reagent.

Lyophilized Platelets are essential in the assessment of von Willebrand Factor (vWF) activity and function through the vW Cofactor Assay. They are reconstituted and used to evaluate the ability of vWF to mediate platelet adhesion and aggregation, critical components of the hemostatic process. When used in conjunction with specific reagents like Ristocetin, Lyophilized Platelets facilitate the von Willebrand Factor Ristocetin Cofactor Assay (vWF), enabling the measurement of vWF's effectiveness.

vW Abnormal Control Plasma is essential for evaluating the functional activity of von Willebrand Factor (vWF) through the vW Cofactor Assay. vW Abnormal Control Plasma, depleted of vWF, is crucial for ensuring the accuracy and reliability of tests by providing a baseline for abnormal vWF activity. The vW Abnormal Control Plasma helps assess the ability of vWF to promote the agglutination of formalin-fixed or fresh platelets in a patient's plasma. A decrease in Ristocetin Cofactor activity, as observed with vW Abnormal Control Plasma, indicates an impairment or dysfunction of vWF.

vW Normal Control Plasma is integral for assessing von Willebrand Factor (vWF) activity and function in various assays. This lyophilized citrated human plasma, confirmed to have a vWF activity level between 80% and 120%, is reconstituted and used to measure the efficacy of vWF in promoting platelet adhesion and aggregation. It serves as a control in the von Willebrand Factor Ristocetin Cofactor Assay (vWF), facilitating the accurate determination of vWF activity in conjunction with reagents such as Ristocetin

vW Normal Reference Plasma is vital for assessing the functional activity of von Willebrand Factor (vWF) in various assays, including the vW Cofactor Assay. This plasma, characterized by a validated vWF activity level within the normal range, serves as a benchmark for evaluating normal vWF activity and ensuring the accuracy and reliability of testing results. vW Normal Reference Plasma facilitates the assessment of vWF's ability to promote the agglutination of platelets, whether formalin-fixed or fresh, in patient plasma. Consistent results within the expected range help confirm the integrity of the assay system, while deviations may indicate potential issues with vWF functionality or assay performance.

TRIS Buffered Saline (TBS) in von Willebrand Factor (vWF) testing plays a crucial role by maintaining a stable pH around 7.4 and adjusting the ionic strength. This facilitates optimal interactions between vWF, Ristocetin, and platelets during the vW Cofactor Assay. TBS creates conditions that mimic physiological shear stress, crucial for assessing vWF's ability to promote platelet adhesion and aggregation. TRIS Buffered Saline (TBS) enhances the reproducibility of assay results across different laboratories, providing reliable insights into vWF function.

AUTOMATION

vW Select™ is intended for use in semi-automated and automated Light Transmission Platelet Aggregometers. This Test System may also be used with other turbidometric or impedance analyzers, and flow cytometers.

QUALITY / QUANTITY

There are no primary standards for vW Select™. The responses to this Test System are concentration dependent. A known normal donor should be tested with each new lot of AggRecetin Reagent, AggRecetin Diluent, Lyophilized Platelets, vW Abnormal Control Plasma, vW Normal Control Plasma, vW Normal Reference Plasma, and TRIS Buffered Saline (TBS). Standards organizations classify Ristocetin Cofactor Activity Assays as semi-quantitative or semi-qualitative.

vW Select™ comes packaged as 2 x 15 mL vials of AggRecetin Reagent, 2 x 2.0 mL vials of AggRecetin Diluent, 2 x 10 mL vials of Lyophilized Platelets, 5 x 0.5 mL vials of vW Abnormal Control Plasma, 5 x 0.5 mL vials of vW Normal Control Plasma, 5 x 0.5 mL vials of vW Normal Reference Plasma, and 3 x 10 mL vials of TRIS Buffered Saline (TBS). The working concentration of Ristocetin is 15 mg / mL.

SPECIMEN TYPE

The test specimen is prepared from sodium citrate anti-coagulated whole blood, and the test sample is Platelet Poor Plasma (PPP). The test blank consists of Lyophilized Platelets and TRIS Buffered Saline (TBS), providing a standardized baseline for the assay.

Ristocetin Reagent may be used with human or animal Platelet Poor Plasma (PPP). Results are based on the extent and rate of platelet agglutination, reflecting the functional activity of vWF.

Reference plasmas, lyophilized formalin-fixed platelets, control plasmas, and diluents are used to perform a Ristocetin Cofactor Activity Test. Test results are determined by interpolation from a standard curve.

TESTING POPULATION

- Human: The prevalence of von Willebrand platelet disorders is global and may vary by race, ethnicity, blood type, and other factors. The incidence is ~2%.
- Anti-Platelet Drugs: The prevalence and incidence are variable. BTK inhibitors and vancomycin are known to decrease RIPA outcomes. A recently developed anti-platelet glycoprotein (GP) Ib monoclonal antibody (moAB) labeled as OP-FI, along with a thoroughly studied anti-GBIb MoAB known as AP-1, completely eliminate platelet agglutination induced by Ristocetin.
- Inherited Platelet Disorders: The prevalence and incidence are variable. Platelets derived from individuals with Bernard-Soulier Syndrome do not agglutinate when exposed to Ristocetin. In contrast to von Willebrand Disease, the levels of von Willebrand Factor activity and von Willebrand antigen remain within normal ranges.
- Animal: The prevalence and incidence are species dependent.

IN VITRO DIAGNOSTIC

vW Select™ contents are in vitro diagnostic reagents intended for Professional Laboratory Use Only. The Test System contents are not intended for injection or ingestion.

INTENDED USER

vW Select™ is intended for Professional Laboratory Use by qualified personnel.

TEST PRINCIPLE

When introduced to a 37°C Platelet Poor Plasma (PPP) test sample, the Ristocetin Reagent induces platelet agglutination by interacting with the glycoprotein Ib (GP Ib) receptor on the platelet surface, which is mediated by von Willebrand Factor (vWF). This initial aggregation, called primary aggregation, is characterized by a shape change in the platelets and is reversible. The primary aggregation may be followed by the release of endogenous ADP from platelet granules, which triggers a secondary, irreversible wave of aggregation. The extent and rate of both primary and secondary aggregation are influenced by the functional activity of vWF in the sample. The Light Transmission Platelet Aggregometer effectively captures these changes by displaying parameters such as the lag phase, shape change, and the rate and extent of aggregation over a predetermined testing period.

CALIBRATORS AND CONTROLS

There are no calibrators or controls required for vW Select™. A known donor sample should be tested with each lot of vW Select™. Responses are concentration dependent.

REAGENT LIMITATIONS

vW Select™ will perform as specified when the Instructions for Use are followed. The Test System must be used prior to the expiration date printed on each vial.

REAGENTS PROVIDED

REF	106730:	2 vials of Ristocetin Reagent (15 mL)
		2 vials of Ristocetin Reagent (2 mL)
		2 vials of Lyophilized Platelets (10 mL)
		5 vials of vW Abnormal Control Plasma (0.5 mL)
		5 vials of vW Normal Control Plasma (0.5 mL)
		5 vials of vW Normal Reference Plasma (0.5 mL)
		3 vials of TRIS Buffered Saline (TBS) (10 mL)

REAGENTS AND MATERIALS REQUIRED BUT NOT PROVIDED

- Purified Water (Distilled, Deionized, Reagent Grade), pH 5.3 – 7.2 for reconstitution


 **NOTE: USING BLOOD BANK SALINE WILL CAUSE ERRONEOUS RESULTS.**


MATERIALS AND ACCESSORIES


- Platelet Aggregometer (Follow the Manufacturer's Instructions for Use)
- Centrifuge
- Electronic Pipette
- Pipette Tips ②
- Aggregometer Test Tubes (Siliconized) ②
- Aggregometer Stir Bars (Plastic Coated) ②
- Plastic Sample Tubes and Caps (for Dilutions) ②
- Mechanical Specimen Rocker


 **NOTE: DISPOSABLE ITEMS SUCH AS TEST TUBES, STIR BARS, SAMPLE TUBES, AND CAPS ARE FOR ONE TIME USE ONLY**


STORAGE AND STABILITY


 AggRecetin Reagent, AggRecetin Diluent, Lyophilized Platelets, vW Abnormal Control Plasma, vW Normal Control Plasma, vW Normal Reference Plasma, and TRIS Buffered Saline (TBS) do not require temperature protection during shipment.

 Upon receipt, store AggRecetin Reagent, AggRecetin Diluent, Lyophilized Platelets, vW Abnormal Control Plasma, vW Normal Control Plasma vW Normal Reference Plasma, and TRIS Buffered Saline (TBS) at 2 – 8° C in their original packaging.


 Resuspended Lyophilized Platelets are stable for 30 days, when stored in its tightly capped, original container at 2 – 8° C.

 Reconstituted AggRecetin Reagent is stable for 7 days, when stored in its tightly capped, original container at 2 – 8° C.


 Reconstituted vW Abnormal Control Plasma, vW Normal Control Plasma, and vW Normal Reference Plasma are stable for 8 hours, when stored in their tightly capped, original container at 2 – 8° C.


 Dilutions containing vW Abnormal Control Plasma, vW Normal Control Plasma, and vW Normal Reference Plasma are stable for 45 minutes at room temperature.


STERILITY

 vW Select™ are not sterile products. Be careful not to contaminate the products when pipetting the reconstituted or aliquoted reagents.


WARNINGS AND PRECAUTIONS


 Wear PPE in accordance with laboratory policies and practices when handling AggRecetin Reagent, AggRecetin Diluent, Lyophilized Platelets, vW Abnormal Control Plasma, vW Normal Control Plasma, vW Normal Reference Plasma, and TRIS Buffered Saline (TBS).


 Follow standard precautions when preparing test specimens and samples.


 Handle AggRecetin Reagent, AggRecetin Diluent, Lyophilized Platelets, vW Abnormal Control Plasma, vW Normal Control Plasma, vW Normal Reference Plasma, and TRIS Buffered Saline (TBS) with care to avoid contamination during use.

 Avoid reagent evaporation by limiting air – liquid exchange surfaces.

 To ensure optimum test results, a known donor control sample should be run consecutively, without interruption.

 To preserve reagent stability, store remaining reagents in their tightly capped, original containers.

 Dispose of post-test materials in accordance with applicable regulations and laboratory policies.

 **NOTE TO USER: ANY SERIOUS INCIDENT THAT OCCURS IN RELATION TO THIS PRODUCT SHALL BE REPORTED TO THE MANUFACTURER AND THE COMPETENT AUTHORITY OF THE MEMBER STATE IN WHICH THE USER AND / OR PATIENT ARE ESTABLISHED.**

INFECTIOUS MATERIAL STATUS

AggRecetin Reagent, AggRecetin Diluent, and Tris Buffered Saline (TBS) do not contain any infectious materials. Lyophilized Platelets, vW Abnormal Control Plasma, vW Normal Control Plasma, and vW Normal Reference Plasma have been screened at the source and confirmed negative for HIV-1 antigen (HIV-1Ag), anti-HIV-1/2 antibodies, Hepatitis B surface antigen (HBsAg), Hepatitis C antibodies, Human T-Lymphotropic Virus Type I and II (anti-HTLV I/II) antibodies, and Syphilis by serological testing. Nevertheless, due to their human origin, all plasma and platelets should be handled as potentially hazardous materials. Test specimens and samples must be considered infectious and should be handled as if capable of transmitting infection. After testing, test specimens and samples must be disposed of in compliance with applicable regulations and laboratory policies.

SPECIAL FACILITIES

vW Select™ does not require the use of special facilities within a laboratory environment.

PREPARATION FOR USE

 **NOTE: vW SELECT CONTENTS MUST BE AT ROOM TEMPERATURE (15–28° C) PRIOR TO RECONSTITUTION. STORED PLATELETS, REAGENTS, AND PLASMAS MUST BE BROUGHT TO ROOM TEMPERATURE PRIOR TO USE.**

vW SELECT

RISTOCETIN COFACTOR ASSAY SYSTEM

PART NUMBER: _____

LOT _____



The vW Select™ Test System is a carefully curated combination of standard components specifically designed for the Ristocetin Cofactor Assay. This unique

configuration is assembled to optimize the assay's performance and reliability when using this test kit. Each element has been selected based on extensive research and testing to ensure the highest standards of accuracy and precision.

The Test System includes tailored instructions that guide users through each step of the assay process, enhancing overall effectiveness. These instructions provide clear guidance to achieve consistent and reproducible results. The vW Select™ Test System comprises specific materials and lot combinations that have undergone rigorous testing to meet elevated performance standards. It is essential to adhere to the specified materials and lot combinations, as changes or substitutions could compromise the assay's integrity and lead to inaccurate results.

RESUSPENSION OF PLATELETS


LYOPHILIZED PLATELETS

CONTENTS: 2 X 10.0 mL Vials

LOT _____ 

To Prepare 10 mL of Lyophilized Platelets

- Add 10 mL of TRIS Buffered Saline (TBS) to the Lyophilized Platelets vial
- Mechanically mix in specimen rocker at room temperature for 30 minutes
- Resuspended Lyophilized Platelets should be kept capped prior to use.

 **NOTE: AFTER REFRIGERATION, MECHANICALLY MIX RESUSPENDED LYOPHILIZED PLATELETS IN ROCKER FOR AT LEAST 30 MINUTES AT ROOM TEMPERATURE TO EQUILIBRATE AND DEGAS THE SUSPENSION PRIOR TO USE.**

While the change in activity of the reconstituted platelet suspension during storage is typically minimal and not significant for standard assays, it is important to monitor the assay system with control materials. This practice ensures that any potential undesirable changes in platelet activity are identified.

Enhanced results can be achieved by storing the reconstituted reagent material on the PAP-8E while stirring in the reagent wells. This method positively impacts the stability and activity of the reagent. The lyophilized platelets can be stored at 37°C with stirring for up to 6 hours. However, any material that has been stored at this temperature for longer than 6 hours should be discarded.

RECONSTITUTION

AGGRECETIN REAGENT

CONTENTS: 2 X 15.0 mL Vials

LOT _____ 

The working concentration of reconstituted AggRecetin Reagent is 15 mg / mL. All final concentrations are based on adding 25 µL of Ristocetin Reagent to a 225 µL Platelet Poor Plasma (PPP) test sample. Bio/Data Corporation has tested the specified lot of AggRecetin Reagent and determined that optimal performance is achieved when a vial is reconstituted with:

_____ mL's of the provided AggRecetin Diluent

- Add the specified volume of AggRecetin Diluent to the AggRecetin Reagent vial
- Invert gently to mix at room temperature
- Reconstituted AggRecetin Reagent should be kept capped prior to use.
- After Reconstitution, allow AggRecetin Reagent stand at room temperature for 30 minutes to rehydrate prior to use.

AGGRECETIN DILUENT

CONTENTS: 2 x 2.0 mL Vials

LOT _____ 

vW ABNORMAL CONTROL PLASMA

CONTENTS: 5 x 0.5 mL Vials

LOT _____ 

- Add 0.5 mL of Purified Water to the vW Abnormal Control Plasma vial
- Invert gently to mix at room temperature
- Reconstituted vW Abnormal Control Plasma should be kept capped prior to use.
- After Reconstitution, allow vW Abnormal Control Plasma stand at room temperature for 20 minutes to rehydrate prior to use.

The accuracy of recovering specific values in the lower range of the assay depends on the quality of the constructed curve. If the assay value falls below the instrument's reportable range, rerun the sample undiluted. The assay value will be half of the reported undiluted value. Users should establish their own specific ranges based on the accepted curve.

vW NORMAL CONTROL PLASMA

CONTENTS: 5 x 0.5 mL Vials

LOT _____ 

- Add 0.5 mL of Purified Water to the vW Normal Control Plasma vial
- Invert gently to mix at room temperature
- Reconstituted vW Normal Control Plasma should be kept capped prior to use.
- After Reconstitution, allow vW Normal Control Plasma stand at room temperature for 20 minutes to rehydrate prior to use.

Different dilutions (1:2 or 1:4) may exhibit varying levels of relative recoverable activity. Each dilution should have its own reference range, which should be used to control the assay system.

vW NORMAL REFERENCE PLASMA

CONTENTS: 5 x 0.5 mL Vials

LOT _____ 

- Add 0.5 mL of Purified Water to the vW Normal Reference Plasma vial
- Invert gently to mix at room temperature
- Reconstituted vW Normal Reference Plasma should be kept capped prior to use.
- After Reconstitution, allow vW Normal Reference Plasma stand at room temperature for 10 minutes to rehydrate
- Invert gently to mix at room temperature
- Allow vW Normal Reference Plasma stand at room temperature for an additional 5 minutes to rehydrate
- Invert gently to mix at room temperature
- The material is now ready for use


TRIS BUFFERED SALINE

CONTENTS: 3 x 10.0 mL Vials

LOT _____ 


PATIENT PREPARATION

Patients should refrain from taking aspirin or using aspirin-containing medications and products, as well as other medications, supplements, or energy drinks known to affect platelet function for 7 – 10 days prior to specimen collection. Ingestion of fatty foods, dairy products, and smoking should be avoided for 12 hours before specimen collection.

 **NOTE: CONSULTATION WITH A PHYSICIAN IS REQUIRED PRIOR TO MAKING ANY MEDICATION CHANGES.**

SPECIMEN COLLECTION

The specimen should be collected with care to avoid stasis, hemolysis, contamination by tissue fluid and exposure to glass. Specimens must be kept at room temperature. Release the tourniquet as soon as blood begins to flow into the collection device.


 **PRACTICE STANDARD PRECAUTIONS THROUGHOUT THE SPECIMEN COLLECTION, SAMPLE PREPARATION, AND ANALYTICAL PROCESSES. DISPOSE OF SHARPS AND BIOHAZARDOUS WASTE IN ACCORDANCE WITH APPLICABLE REGULATIONS AND LABORATORY POLICIES.**

Evacuated Specimen Collection Technique

- Use a 21g or 23g winged needle collection set for specimen collection
- Draw blood into plastic evacuated specimen collection tubes containing 3.2% (0.11 M) sodium citrate anti-coagulant
- Gently mix the specimen collection tube 4 - 5 times by inversion
- Write collection time on the specimen label
- Maintain specimen collection tubes at room temperature
- Remix specimen collection tubes prior to centrifugation

Syringe Collection Technique

- Use a 21g or 23g winged needle collection set for the venipuncture
- Draw 9.0 mL of blood into a plastic syringe, avoiding excess suction
- Clamp the winged needle tubing and disconnect the syringe
- Immediately and gently dispense the blood specimen into a plastic (polypropylene) tube containing 1.0 mL of 0.11 M sodium citrate anti-coagulant. The blood to anticoagulant ratio is 9 parts blood to 1 part anti-coagulant
- Cap the plastic tube
- Gently mix the specimen collection tube 4 - 5 times by inversion
- Write collection time on the specimen label
- Maintain specimen collection tubes at room temperature
- Remix specimen collection tubes prior to centrifugation

 **NOTE: WHEN THE PATIENT'S HEMATOCRIT IS LESS THAN 30% OR GREATER THAN 55%, THE BLOOD TO ANTI-COAGULANT RATIO MUST BE ADJUSTED. BLUE TOP EVACUATED SPECIMEN COLLECTION TUBES MUST CONTAIN 3.2% (0.11 M) SODIUM CITRATE ANTICOAGULANT, WHICH IS THE RECOMMENDED CONCENTRATION FOR PLATELET FUNCTION STUDIES.**

SAMPLE PREPARATION

Platelet Rich Plasma (PRP)

- Centrifuge the anti-coagulated blood at 150 x g for 10 minutes at room temperature
- Examine the plasma layer for red cells
- If red cells are present, re-centrifuge for an additional 5 minutes
- Use a Pipette to transfer the PRP to a plastic container labeled PRP
- Remove the PRP from a point just below the middle of the PRP volume for consistent platelet count (*THE TOP OF THE VOLUME HAS A LOWER PLATELET COUNT AND THE BOTTOM IS MORE CONCENTRATED*)
- Cap the container
- Allow the container to stand at room temperature

Platelet Poor Plasma (PPP)

- Centrifuge the remaining PRP blood specimen at 2500 x g for 20 minutes
- Use a Pipette to transfer the PPP to a plastic container labeled PPP
- Cap the container
- Allow the container to stand at room temperature

ASSAY PROCEDURE

Ristocetin Cofactor Assay Three Point Standard Curve Procedure



NOTE: THIS IS A GENERAL PROCEDURE. FOLLOW THE INSTRUCTIONS FOR USE PROVIDED BY THE MANUFACTURER OF THE AGGREGOMETER IN USE.

Prepare a Blank for Each Patient



NOTE: IN RISTOCETIN COFACTOR ASSAY TESTING, THE SAME BLANK MAY BE USED FOR EACH TEST WELL IN BOTH THE RISTOCETIN COFACTOR ASSAY THREE POINT STANDARD CURVE AND THE RISTOCETIN COFACTOR ASSAY THREE POINT PATIENT TESTING, PROVIDED THAT ALL TESTING IS COMPLETED WITHIN FOUR HOURS OF PREPARING THE BLANK.

* _____ Blank TBS Adjustment Factor

To achieve optimal performance, a precise ratio of Lyophilized Platelets and TRIS Buffered Saline (TBS) will be used to prepare the blank. This ratio is based on a fixed volume of 175 μ L of Lyophilized Platelets and a variable volume of TRIS Buffered Saline (TBS). The volume ratios may be adjusted to accommodate the requirements of the aggregometer.

The Blank TBS Adjustment Factor*, provided above, is multiplied by the fixed volume of 175 μ L of Lyophilized Platelets to determine the appropriate amount of TRIS Buffered Saline (TBS) needed for preparing the blank.

For example, if the TBS Adjustment Factor is 0.9, the amount of TBS to be added to the fixed volume of 175 μ L of Lyophilized Platelets will be 158 μ L. The calculation is $0.9 \times 175 \mu\text{L} = 157.5 \mu\text{L}$, which is rounded up to 158 μ L.

$$\text{_____ TBS Adjustment Factor} \times 175 \mu\text{L of Lyophilized Platelets} \\ = \text{_____ } \mu\text{L of TBS, round up to _____ } \mu\text{L}$$

- Label a test tube with the letter "M", test well #, and patient ID to identify the Blank Mixture
- Pipette 175 μ L of Lyophilized Platelets and Amount Calculated Above of Tris Buffered Saline (TBS) into the test tube (*DO NOT ADD A STIR BAR*)
- Invert gently to mix
- Label a test tube with the letter "B" to identify the Blank
- Add a Stir Bar to test tube labeled "B"
- Pipette 250 μ L of the test tube labeled "M" (Blank mixture) to the test tube labeled "B"
- Discard the test tube labeled "M"
- Seal test tube labeled "B" with Parafilm® or similar material
- Place test tube labeled "B" aside for later use



NOTE: THE PREPARED BLANK IS STABLE FOR 4 HOURS AT ROOM TEMPERATURE

Prepare Standard Curve Dilutions



NOTE: DILUTIONS SHOULD BE PREPARED INDEPENDENTLY AND NOT AS SERIAL DILUTIONS.

Use the Following Dilution Procedure for a 100%, 50%, and 25% Curve:

- Label three new test tubes: 100%, 50%, and 25%
- In the test tube labeled 100%: Pipette 200 μ L of vW Normal Reference Plasma (NRP) and 200 μ L of Tris Buffered Saline (TBS)
- Invert gently to mix
- In the test tube labeled 50%: Pipette 100 μ L of vW Normal Reference Plasma (NRP) and 300 μ L of Tris Buffered Saline (TBS)
- Invert gently to mix
- In the test tube labeled 25%: Pipette 50 μ L of vW Normal Reference Plasma (NRP) and 350 μ L of Tris Buffered Saline (TBS)
- Invert gently to mix
- Let each dilution stand for 10 minutes
- Place all three vW Normal Reference Plasma dilutions aside for later use

Table 1: Standard Curve Dilutions for Ristocetin Cofactor Testing

DILUTIONS	VOLUME OF NRP	VOLUME OF TBS
100% CURVE (1:2)	200	200
50% CURVE (1:4)	100	300
25% CURVE (1:8)	50	350
12.5% CURVE (1:16) OPTIONAL	25	375



NOTE: DILUTIONS ARE STABLE FOR 40 MINUTES AT ROOM TEMPERATURE

Prepare Samples

- Label three new test tubes with 100% Well #1, 50% Well # 2, and 25% Well # 3
- Place the labeled test tubes into the corresponding well # 1 - 3 of the Stirred Sample Incubation Wells
- Add a Stir Bar to Each Test Tube
- Pipette 200 μ L of Lyophilized Platelets into each Test tube in the stirred sample incubation wells (*MAKE SURE THERE ARE NO BUBBLES*)
- Pipette 25 μ L of Ristocetin Reagent into each test tube in the stirred sample incubation wells (*MAKE SURE THERE ARE NO BUBBLES*)
- Select the on-screen timer for each stirred sample incubation wells and the warming countdown will start
- The samples will incubate at 37°C for the pre-set time

Set the 100% Baseline (Blank)

- Locate the previously prepared Blank test tube labeled "B"
- Place the test tube into test well # 1
- Select BLANK to activate the test well
- The BLANK button will change to START
- Repeat the above process for test wells # 2 and # 3

Begin Testing

- Once the countdown timer reaches 0:00, press the timer button to stop each stirred sample incubation wells
- Transfer the test tube in the stirred sample incubation well #1 to test well # 1
- Repeat the step above for each test well, making sure all test tubes remain with their corresponding well #'s during transfer
- Close the pipette guides
- Select START for test well # 1
- Pipette 25 μ L of the previously prepared dilution labeled 100% into test well # 1 (*DO NOT ALLOW REAGENT TO RUN DOWN THE WALL OF THE TEST TUBE OR PERMIT THE PIPETTE TIP TO BREAK THE SURFACE OF THE SAMPLE*)
- Select INJECT for test well # 1
- Select START for test well # 2
- Pipette 25 μ L of the previously prepared dilution labeled 50% into test well # 2 (*DO NOT ALLOW REAGENT TO RUN DOWN THE WALL OF THE TEST TUBE OR PERMIT THE PIPETTE TIP TO BREAK THE SURFACE OF THE SAMPLE*)
- Select INJECT for test well # 2
- Select START for test well# 3
- Pipette 25 μ L of the previously prepared dilution labeled 25% into test well# 3 (*DO NOT ALLOW REAGENT TO RUN DOWN THE WALL OF THE TEST TUBE OR PERMIT THE PIPETTE TIP TO BREAK THE SURFACE OF THE SAMPLE*)
- Select INJECT for test well # 3
- The test will now run for the pre-set time (*OTHER MANUFACTURER'S TEST PROCEDURES MAY SPECIFY DIFFERENT TIMES OR VOLUMES*)

Ristocetin Cofactor Assay Three Point Patient Test Procedure



NOTE: THIS IS A GENERAL PROCEDURE. FOLLOW THE INSTRUCTIONS FOR USE PROVIDED BY THE MANUFACTURER OF THE AGGREGOMETER IN USE.

Prepare a Blank for Each Patient



NOTE: IN RISTOCETIN COFACTOR ASSAY TESTING, THE SAME BLANK MAY BE USED FOR EACH TEST WELL IN BOTH THE RISTOCETIN COFACTOR ASSAY THREE POINT CURVE AND THE RISTOCETIN COFACTOR ASSAY THREE POINT PATIENT TESTING, PROVIDED THAT ALL TESTING IS COMPLETED WITHIN FOUR HOURS OF PREPARING THE BLANK.

- Use the Blank previously prepared for the Ristocetin Cofactor Assay Three Point Curve
- If the Blank was prepared more than 4 hours ago, it is no longer stable, and a new Blank must be prepared
- To prepare a new Blank, follow the blank preparation instructions located in this IFU under Ristocetin Cofactor Assay Three Point Standard Curve Procedure

Prepare Dilutions

- Label one to eight new test tubes with each patient's ID and well #
- Pipette 100 μ L of the patient sample into the corresponding labeled test tube for each patient being tested
- Pipette 100 μ L of TRIS Buffered Saline (TBS) into each test tube
- Place dilutions aside for later use

Prepare Samples

- Label one to eight new test tubes with each patient's ID and well #
- Place the labeled test tubes into the corresponding wells # 1 - 8 of the stirred sample incubation wells
- Add a Stir Bar to each test tube
- Pipette 200µL of Lyophilized Platelets into each test tube in the stirred sample incubation wells (MAKE SURE THERE ARE NO BUBBLES)
- Pipette 25µL of Ristocetin into each test tube in the stirred sample incubation wells (MAKE SURE THERE ARE NO BUBBLES)
- Select the on-screen timer for each stirred sample incubation wells and the warming countdown will start
- The samples will incubate at 37°C for the pre-set time

Set the 100% Baseline (Blank)

- Locate the previously prepared Blank test tube labeled "B"
- Place the test tube into test well # 1
- Select BLANK to activate the test well
- The BLANK button will change to START
- Repeat the above process for each test well being used for testing

Begin Testing

- Once the countdown timer reaches 0:00, press the timer button to stop each stirred sample incubation wells
- Transfer the test tube from the stirred sample incubation wells # 1 to test well # 1
- Repeat the step above for each test well, making sure all test tubes remain with their corresponding cell #'s during transfer
- Close the pipette guides
- Select START for test well # 1
- Pipette 25µL of the previously prepared patient dilution directly into the test tube in test well # 1
- VERIFY THE CORRECT PATIENT DILUTIONS ARE BEING TRANSFERRED TO THE CORRECT PATIENT TEST WELL (DO NOT ALLOW REAGENT TO RUN DOWN THE WALL OF THE TEST TUBE OR PERMIT THE PIPETTE TIP TO BREAK THE SURFACE OF THE SAMPLE)
- Select INJECT for test well # 1
- Repeat the above procedure for each test well being used for testing
- The test will now run for the pre-set time (OTHER MANUFACTURER'S TEST PROCEDURES MAY SPECIFY DIFFERENT TIMES OR VOLUMES)

NOTE: USE A KNOWN DONOR AS A CONTROL SAMPLE. EACH LABORATORY SHOULD ESTABLISH AND VALIDATE ITS OWN TEST PROTOCOL AND VERIFY THE RESULTING PERFORMANCE OF ITS TEST SYSTEM (REAGENTS, INSTRUMENT, AND TEST PROTOCOL).

QUALITY CONTROL

For platelet aggregation studies, a known donor should be tested in the same manner as the patient to ensure test system performance and consistency. A new control should be included with each test series, and preferably with each new reagent lot or after instrument maintenance. Each laboratory must define its acceptable ranges for its patient population and verify the expected performance of the test system.

A von Willebrand Factor deficient plasma is included as an abnormal control and should be assayed as a test plasma with an expected result of $\leq 45\%$ activity. This control ensures that the assay system is specific for the von Willebrand Factor and that agglutination will not be influenced by other normal plasma proteins. Additionally, it is suggested to run vW Normal Reference Plasma and vW Abnormal Control Plasma to validate standard curves.

RESULTS

The vW Ristocetin Cofactor Assay utilizes a three-point standard curve to determine the percent activity of vWF in a patient's plasma. As illustrated in Figure 1, the standard curve is generated using three known dilutions of vWF, typically representing 100%, 50%, and 25% activity. The assay measures the response at these concentrations and plots the values to establish a reference curve.

When the patient's sample is analyzed, its response is interpolated against the standard curve to determine vWF activity as a percentage of normal plasma (100% activity). Clinically, vWF activity levels within the normal range (45–150%) indicate typical function, while values below 45% may suggest von Willebrand Disease (vWD). Conversely, elevated vWF activity may be associated with conditions or disorders.

Figure 2 illustrates key phases of platelet aggregation over time, providing insights into platelet function. The initial downward deflection observed in each curve represents shape change, where platelets respond to the agonist by undergoing morphological alterations without immediate aggregation. This is followed by the primary aggregation phase, characterized by an upward trend as platelets begin to clump together, increasing light transmission. The maximum aggregation point varies across the curves, reflecting differences in platelet responsiveness. Certain tracings, such as those in channels 6, 7, and 8, exhibit a decline after reaching maximum aggregation, indicating a reversible aggregation pattern, whereas others sustain aggregation, signifying a strong response. Variations in lag time before aggregation initiation highlight differences in platelet activation kinetics. These characteristics help assess platelet function and identify potential abnormalities in aggregation response.

The von Willebrand Factor (vWF) Assay tracing shown in Figure 3, illustrates an abnormal vWF activity response, with results measuring below 45% activity. This reduced activity suggests a potential von Willebrand Disease (vWD). The tracing demonstrates a weaker-than-expected response compared to the standard curve, indicating impaired vWF function or decreased factor levels. Such abnormal results may be associated with vWD Type 1 or Type 2 variants, acquired vWF deficiency, or other coagulation disorders. These findings underscore the need for further diagnostic testing to accurately characterize the nature and severity of the deficiency.

It is essential to interpret these aggregation results within the broader context of the patient's clinical condition. A definitive diagnosis should only be made after further testing and comprehensive evaluation. The figures include spike marks that indicate the precise points of reagent addition, providing clear reference points for understanding the timing of reagent introduction and its immediate effects on the aggregation process.

FIGURE 1: VW RISTOCETIN COFACTOR ASSAY THREE POINT PERCENT ACTIVITY




FIGURE 2: NORMAL VW RISTOCETIN COFACTOR ASSAY THREE POINT PATIENT TESTING





FIGURE 3: ABNORMAL VW RISTOCETIN COFACTOR ASSAY THREE POINT PATIENT TESTING



EXPECTED VALUES

Each laboratory should establish expected ranges for each reagent at various concentrations used to induce aggregation.

A result of less than 40% von Willebrand Factor is considered abnormal and suggestive of von Willebrand Syndrome. However, other properties of the von

Willebrand molecule must be considered for diagnosis of the variant forms of von Willebrand Syndrome. Since reference ranges for von Willebrand Factor are dependent on blood type, each laboratory should establish blood type specific reference ranges for its patient population.

vW Abnormal Control Plasma will yield von Willebrand Factor Assay results of $\leq 45\%$. The ability to generate a quantitative value in this range is dependent upon the sensitivity of the assay system in use.

LIMITATIONS

In Light Transmission Aggregometry, the presence of red blood cells in the PRP will cause the observed aggregation to be reduced. The presence of platelets in the PPP will cause final aggregation to be increased. Spurious results may occur if the PRP platelet count is less than 75,000 platelets / cumm. PRP platelet counts can only be performed using the hemocytometer method. Compromised samples must be rejected. If the results are abnormal, the test should be repeated on another occasion. Each laboratory must establish reference ranges tailored to the population it serves, and the specific reagent concentrations used.

ANALYTICAL PERFORMANCE

Platelet aggregation, induced by commonly used reagents like Ristocetin Reagent, is a non-linear test system. Responses are based on the difference between the patient's Platelet Rich Plasma and Platelet Poor Plasma light transmission and therefore, results are unique to that patient. Certain parameters are more prone to non-linearity than others. These include lag phase, primary slope, secondary slope, biphasic response and disaggregation. The non-linearity is caused by many factors such as the reaction chemistry and instrumentation. Platelet aggregation displays the response rate or activity and does not quantify the reactants or their concentrations.

In platelet aggregation, accuracy is a relative parameter and is dependent on the test system. The limitations of platelet aggregation make it difficult to provide typical precision or reproducibility ranges.

The variability in linearity, precision and reproducibility of results in Ristocetin Reagent-based test systems is acknowledged by multiple standards organizations. The commonly accepted CV is $\pm 15\%$.

Test to Test Reproducibility:	less than $\pm 7.5\%$
Instrument to Instrument Reproducibility:	less than $\pm 15.0\%$
Reagent Lot to Lot Variability:	less than $\pm 10.5\%$
Laboratory to Laboratory (System to System)	less than $\pm 12.5\%$

REFERENCES

- Allain JP, Cooper HA, Wagner RH, Brinkhous KM. Platelets fixed with paraformaldehyde: a new reagent for assay of von Willebrand factor and platelet aggregating factor. *J Lab Clin Med.* 1975 Feb;85(2):318-28.
- Angiolillo DJ, Ueno M, Goto S. Basic principles of platelet biology and clinical implications. *Circ J.* 2010 Apr;74(4):597-607.
- Born GV, Cross MJ. The Aggregation of Blood Platelets. *J Physiol.* 1963 Aug; 168(1):178-95.
- Brinkhous KM, Graham JE, Cooper HA, Allain JP, Wagner RH. Assay of von Willebrand factor in von Willebrand's disease and hemophilia: use of a macroscopic platelet aggregation test. *Thromb Res.* 1975 Mar;6(3):267-72.
- Brinkhous KM, Read MS. Preservation of platelet receptors for platelet aggregating factor/von Willebrand factor by air drying, freezing, or lyophilization: new stable platelet preparations for von Willebrand factor assays. *Thromb Res.* 1978 Oct;13(4):591-7.
- Bye A, Lewis Y, O'Grady J. Effect of a single oral dose of aspirin on the platelet aggregation response to arachidonic acid. *Br J Clin Pharmacol.* 1979 Mar; 7(3):283-6.
- Cattaneo M, Cerletti C, Harrison P, Hayward CP, Kenny D, Nugent D, Nurden P, Rao AK, Schmaier AH, Watson SP, Lussana F, Pugliano MT, Michelson AD. Recommendations for the Standardization of Light Transmission Aggregometry: A Consensus of the Working Party from the Platelet Physiology Subcommittee of SSC/ISTH. *J Thromb Haemost.* 2013 Apr 10.
- CLSI. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. CLSI document H18-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
- CLSI. Protection of Laboratory Workers from Occupationally Acquired Infections, Approved Guideline - Fourth Edition. CLSI document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
- CLSI. Platelet Function Testing by Aggregometry, Approved Guideline - Fourth Edition. CLSI document H58-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- CLSI. Collection, Transport and Processing for Plasma Based Coagulation Assays and Molecular Hemostasis Assays, Approved Guideline - Fifth Edition. CLSI document H21-A5. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
- CLSI. Clinical Laboratory Safety, Approved Guideline - Third Edition. CLSI document GP17-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- Day HJ, Holmsen H. Laboratory tests of platelet function. *Ann Clin Lab Sci* (1971). 1972 Jan-Feb; 2(1):63-74.
- Day HJ, Rao AK. Evaluation of platelet function. *Semin Hematol.* 1986 Apr;23(2):89-101.
- Eichelberger, JW. Kinetic (Slope) Measurement of Platelet Aggregation. Bio/Data Corporation, Horsham, PA; 1984.
- Favaloro EJ, Gosselin RC, Pasalic L, Lippi G. Post-analytical issues in hemostasis and thrombosis testing: An update. In EJF, RCG, editors, Hemostasis and Thrombosis: Methods and Protocols. 2nd ed. New York: Humana Press. 2023. p. 787-811. (Methods in Molecular Biology).
- Federici AB, Lee CA, Berntorp EE, Lillicrap D, Montgomery RR. Von Willebrand Disease: Basic and Clinical Aspects. 2011.
- Garner JS. Guideline for isolation precautions in hospitals. The Hospital Infection Control Practices Advisory Committee. *Infect Control Hosp Epidemiol.* 1996 Jan;17(1):53-80.
- Grainick HR, Sultan Y, Collier BS. Von Willebrand's disease: combined qualitative and quantitative abnormalities. *N Engl J Med.* 1977 May 5;296(18):1024-30.
- Harmening, D. M. Clinical Hematology and Fundamentals of Hemostasis. Fifth Edition. F. A. Davis Company. 2009.
- Hoffbrand, A. V., Moss, P. A. H., & Pettit, J. E. Hoffbrand's Essential Haematology. Seventh Edition. John Wiley & Sons Ltd. 2016.
- Howard MA, Firkin BG. Ristocetin--a new tool in the investigation of platelet aggregation. *Thromb Diath Haemorrh.* 1971 Oct 31; 26(2): 362-9.
- Israels SJ, El-Ekiaby M, Quiroga T, Mezzano D. Inherited disorders of platelet function and challenges to diagnosis of mucocutaneous bleeding. *Haemophilia.* 2010 Jul;16 Suppl 5:152-9.
- Kabayashi J, Shinoki N, Nakamura T, Ariyoshi H, Kawasaki T, Sakon M, Monden M. Prevalence of impaired responsiveness to epinephrine in platelets among Japanese. *Thromb Res.* 1996 Jan 1;81(1):85-90.
- Kaushansky K, Lichtman MA, Prchal JT, Levi MM, Press OW, Burns LJ, Caligiuri M. eds. Williams Hematology, 9e. McGraw-Hill Education. 2015.
- Keohane, E. M., Smith, L. J., Walenga, J. M., & Block, D. R. Rodak's Hematology: Clinical Principles and Applications. Fifth Edition. Saunders, an imprint of Elsevier Inc. 2016.
- Levine PH. The effect of thrombocytopenia on the determination of platelet aggregation. *Am J Clin Pathol.* 1976 Jan;65(1):79-82
- Linnemann B, Schwonberg J, Mani H, Prochnow S, Lindhoff-Last E. Standardization of light transmittance aggregometry for monitoring antiplatelet therapy: an adjustment for platelet count is not necessary. *J Thromb Haemost.* 2008 Apr;6(4):677-83.
- Marcus AJ, Coleman RW, Hirsh J, Ivarer VJ, Salzman EW. Hemostasis and thrombosis: Basic Principles and Clinical Practice. Vol. 472. Philadelphia: JB Lippincott Company; 1982.
- Michelson, AD. Platelets. Third Edition. Amsterdam: Academic Press; 2013.
- Miller CH, Graham JB, Goldin LR, Elston RC. Genetics of classic von Willebrand's disease. I. Phenotypic variation within families. *Blood.* 1979 Jul;54(1):117-36.
- Mills DC, Robb IA, Roberts GC. The release of nucleotides, 5-hydroxytryptamine and enzymes from human blood platelets during aggregation. *J Physiol.* 1968 Apr;195(3):715-29.
- Moncada S, Vane JR. Arachidonic acid metabolites and the interactions between platelets and blood-vessel walls. *N Engl J Med.* 1979 May 17;300(20):1142-7.
- NCCLS. Assays of von Willebrand Factor Antigen and Ristocetin Cofactor Activity; Approved Guideline. NCCLS document H51-A. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2002.
- Nilsson, I. M. and Holmberg, L.: von Willebrand's Disease Today. Clin. Hematol., 8:276, 1979.
- O'Donnell CJ, Larson MG, Feng D, Sutherland PA, Lindpaintner K, Myers RH, D'Agostino RA, Levy D, Tofer GH; Framingham Heart Study. Genetic and environmental contributions to platelet aggregation: the Framingham heart study. *Circulation.* 2001 Jun 26;103(25):3051-6.
- Olson JD, Brockway WJ, Fass DN, Magnuson MA, Bowie EJ. Evaluation of ristocetin-Willebrand factor assay and ristocetin-induced platelet aggregation. *Am J Clin Pathol.* 1975 Feb;63(2):210-8.
- Owen CA Jr, Bowie EJW, Thompson JH Jr. The Diagnosis of Bleeding Disorders. 2nd ed. Little, Brown, and Company; 1975.
- Palma-Barqueros V, Revilla N, Sánchez A, Zamora Cánovas A, Rodríguez-Alén A, Marín-Quilez A, González-Porras JR, Vicente V, Lozano ML, Bastida JM, Rivera J. Inherited Platelet Disorders: An Updated Overview. *Int J Mol Sci.* 2021 Apr 26;22(9):4521.
- Ramsey R, Evatt BL. Rapid assay for von Willebrand factor activity using formalin-fixed platelets and microtitration technic. *Am J Clin Pathol.* 1979 Dec;72(6):996-9.
- Siegel JD, Rhinehart E, Jackson M, Chiarello L; Health Care Infection Control Practices Advisory Committee. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Health Care Settings. *Am J Infect Control.* 2007 Dec;35(10 Suppl 2):S65-164.
- The Hospital Infection Control Practices Advisory Committee, Centers for disease Control and Prevention, Public Health Service, U.S. Department of Health and Human Services. Guideline for isolation precautions in hospitals Part II. Recommendations for isolation precautions in hospitals. *American Journal of Infection Control.* 1996; Vol 24, Issue 1: 32-52.
- Triplet DA, et al. Platelet function: laboratory evaluation and clinical application. Chicago, IL: American Society for Clinical Pathology 1978.
- Weiss HJ. Aspirin and Platelets in Drugs and Hematologic Reactions. New York, NY: Dimittov and Nodine, eds. Grune and Stratton. 1974.
- White, M.M., and Jennings, L.K. Platelet Protocols: Research and Clinical Laboratory Procedures, Academic Press, Inc.; 1999.
- Williams WJ, Beutler E, Erslev AJ, Rundles RW. Hematology. New York, NY: McGraw-Hill. 1977.
- Zimmerman TS, Abildgaard CF, Meyer D. The factor VIII abnormality in severe

- von Willebrand's disease. N Engl J Med. 1979 Dec 13;301(24):1307-10.
- Zuzel M, Nilsson IM, Aberg M. A method for measuring plasma ristocetin cofactor activity. Normal distribution and stability during storage. Thromb Res. 1978 May;12(5):745-54.

REVISION HISTORY

Document No: 106826 Revision: AA, August 2025

- Modified Testing Instructions
- Implemented IVDR Regulatory Requirements
- Reformatted and Reconfigured to Enhance Operator Use

For a complete product catalog, please visit our website at www.biodatacorp.com or contact our Customer Service Department.

THE BIO/DATA CORPORATION PRODUCT LINE INCLUDES GENERAL PURPOSE, PROFESSIONAL LABORATORY USE REAGENTS INTENDED TO INDUCE AND REPORT PLATELET FUNCTION ACTIVITY AND RESPONSES. THIS PRODUCT IS WARRANTED TO PERFORM AS DESCRIBED IN ITS LABELING INCLUDING THE INSTRUCTIONS FOR USE. BIO/DATA CORPORATION MAKES NO CLAIM OR WARRANTY, EXPRESSED OR IMPLIED, OF THE CAPABILITY, FITNESS, OR MERCHANTABILITY FOR ANY OTHER PURPOSE. IN NO EVENT SHALL BIO/DATA CORPORATION BE LIABLE FOR ANY CONSEQUENTIAL DAMAGES ARISING OUT OF AFORESAID EXPRESSED WARRANTY.

SYMBOLS



Bio-Hazardous



Catalog Number



Caution



CE Marked & Registered Product



Consult Instructions For Use



European Union Representative



In Vitro Diagnostic Device



Manufacturer



Must Read



Non-Sterile



Single Use Only



Temperature Limitations



United Kingdom Marked & Registered Product



United Kingdom Representative