

## Stir Bars & Test Tubes In One!

**Now  
Conveniently  
Packaged  
Together**



Each Micro Test Tube Contains One Pre-Loaded Stir Bar.



- Eliminates Operator Error
- Saves Preparation Time
- Easy To Handle
- Prevents Contamination
- Economically Priced

MagneTube™ combines two essential products needed for platelet aggregation: siliconized test tubes and plastic coated stir bars. This allows the operator more time to focus on testing rather than on preparation.

Each 7.25 x 55mm siliconized micro test tube contains one pre-loaded micro stir bar for use with the PAP-8E & the PAP-4 series aggregometers using micro-volume adapters.

**NEW PRODUCT**

**Catalog No.  
107146**

**50 Sets per Box**

### BIO/DATA CORPORATION

1-800-257-3282 USA / 215-441-4000 Worldwide

Website: [www.biodatacorp.com](http://www.biodatacorp.com)

Email: [customer.service@biodatacorp.com](mailto:customer.service@biodatacorp.com)

An ISO 9001:2008, 13485:2003 Registered Company

### ORDER TODAY

**Authorized Representative:**  
EMRGO EUROPE, Molenstraat 15, 2513 BH,  
The Hague, The Netherlands

