

Hematology Reagents

Catalog No. 100989, 100994, 100991

Formulated to Enhance Platelet Aggregation & Coagulation Studies



**Packaged for
Use in Procedures
Requiring Different
Volumes**

Calcium Chloride Solution, 0.025 M

Catalog No. 100989
Net Contents: 473 mL

Calcium Chloride Solution is for use in conjunction with other IVD reagents in the performance of coagulation and other laboratory tests, such as PTs, APTTs, factor assays and plasma recalcification (plasma clotting) times.

Sodium Citrate Solution, 0.11 M

Catalog No. 100994
Net Contents: 237 mL

Sodium Citrate Solution is routinely used to anti-coagulate whole blood specimens from human and non-human species for hemostasis testing. It is formulated at 0.11M as recommended by the International Committee of Thrombosis and Hemostasis, Subcommittee on Platelets.

Tris Buffered Saline

Catalog No. 100991
Net Contents: 473 mL

Tris Buffered Saline is routinely used for the reconstitution and dilution of reagents and samples tested for platelet aggregation, Ristocetin CoFactor Activity Assays or other laboratory studies.



155 Gibraltar Road
Horsham, PA 19044 USA

Worldwide: +1 215-441-4000
USA: 1-800-257-3282
Fax Worldwide: +1 215-443-8820
customer.service@biodatacorp.com

REF
100989
100994
100991



An ISO 13485 Registered Company
www.biodatacorp.com



mdi Europa GmbH
Langenhagener Str. 71
D-30855 Langenhagen
GERMANY



Alpha Laboratories
40 Parham Drive Eastleigh
SO50 4NU Hampshire
UNITED KINGDOM

