

# SpeciSafe®

## The Smarter Way to Pack Biological Samples

All-in-one secondary packaging for effortless UN3373 compliance.

Make sample packaging faster, cleaner, and more reliable with SpeciSafe® – the ready to use system that combines absorbency, protection, and tube separation in one robust pack.

### Key Benefits

- ✓ **All in One Design** – Bonded ultra-absorbent layer, built in tube separation, and rigid clam shell case.
- ✓ **Leak-Proof Protection** – Blister type casing fully encloses the primary vessel for secure transport.
- ✓ **No More Wrapping Vials** – Integrated gel block absorbent layer removes the need for separate absorbent materials.
- ✓ **Effortless Compliance** – Supports UN3373 Biological Substance, Category B transport when used with a 95kPa primary container.
- ✓ **Versatile Range** – Multiple formats available to fit the most used 95kPa sample vessels.

### Who Uses SpeciSafe®?

- › Hospitals
- › GP Practices
- › Diagnostic Laboratories
- › Research Facilities
- › Community Testing Services




Packing Instructions  
1. Place primary  
2. Firmly close  
3. Seal secure




# How SpecSafe® Works

## 1 Pack



Push the primary container into pre-moulded shape so they are held in position.

## 2 Close



Close SpecSafe® and work around the edge pressing firmly to seal all around.

## 3 Send



Place SpecSafe® into an appropriate outer packaging.

## Why Choose SpecSafe®?

- › Reduces packaging time
- › Minimises material waste
- › Lowers risk of packaging errors
- › Improves workflow efficiency
- › Enhances sample protection during transit

## Simplify Your Sample Packaging Today



SpecSafe® gives healthcare and laboratory teams a smarter, more efficient way to prepare patient samples, without compromising safety or compliance. Want to learn more or request samples?

Visit: [www.alphalabs.co.uk/specisafe](http://www.alphalabs.co.uk/specisafe) or scan the QR code.

Product Code	Product Description	Pack Size
SS0500	1-place SpecSafe® for Venous Blood Tubes	500
SS0401	2-place SpecSafe® for Venous Blood Tubes	300
SS0701	5-place SpecSafe® for Venous Blood Tubes	150
SS1700	1-place SpecSafe® for Micro/Serum Tubes	450
SS0301	6-place SpecSafe® for Micro/Serum Tubes	300
SS0600	1-place SpecSafe® for Universal Containers	300
SS1400	3-place SpecSafe® for Universal Containers	160

For more information or to purchase,  
please contact Bio/Data at:  
[customer.service@biodatacorp.com](mailto:customer.service@biodatacorp.com)  
or call 1-800-257-3282 (USA)  
+1 215-441-4000 (Worldwide)