

MagneTube™

Stir Bars & Test Tubes In One

Catalog No. 107146

**Conveniently
Packaged
Together**

Each Micro Test Tube Contains One Pre-Loaded Stir Bar

MagneTube™ combines two essential products needed for platelet aggregation: siliconized test tubes and plastic coated stir bars. This allows the operator more time to focus on testing rather than on preparation.



Each 7.25 x 55mm siliconized micro test tube contains one pre-loaded micro stir bar for use with the PAP-8E & the PAP-4 Series aggregometers using micro-volume adapters.



Net Contents: 50 Sets per Box

- **Eliminates Operator Error**
- **Saves Preparation Time**
- **Easy To Handle**
- **Prevents Contamination**
- **Economically Priced**

An ISO 13485:2003 Registered Company

Worldwide: +1 215-441-4000
USA: 1-800-257-3282
Fax Worldwide: +1 215-443-8820
customer.service@biodatacorp.com