

The Gold Standard In Routine Platelet Function Testing

Lyophilized Platelets are a standardized and fixed suspension of human platelets routinely used as a component of a Ristocetin CoFactor Assay Activity Test.



AVAILABLE IN TWO SIZES

Net Contents	Catalog No. 101595	Catalog No. 101258
Lyophilized Platelets	3 x 4.0 mL	1 x 10.0 mL
Tris Buffered Saline	3 x 4.0 mL	1 x 10.0 mL

Platelets You Can Count On

- Standardized Suspension of Formalin – Fixed Human Platelets
- Displays Presence of vWF Inhibition in Mixed Patient and Pooled Plasma Samples
- Retains Gp1b Receptor Function
- Used in Ristocetin CoFactor Activity Assay
- Suitable for Veterinary Use

Lyophilized Platelets are a standardized and fixed platelet suspension derived from human blood. The platelets are lyophilized to ensure stability during long term storage. Pre-measured Tris Buffered Saline Diluent is supplied with the platelets. Upon reconstitution, the platelet suspension will have a count of approximately 200,000/ μ L (using hemocytometer methodology).


 155 Gibraltar Road
 Horsham, PA 19044 USA

Worldwide: +1 215-441-4000
 USA: 1-800-257-3282
 Fax Worldwide: +1 215-443-8820
 customer.service@biodatacorp.com

REF
 101595
 101258


An ISO 13485 Registered Company
www.biodatacorp.com


 mdi Europa GmbH
 Langenhagener Str. 71
 D-30855 Langenhagen
 GERMANY

 Alpha Laboratories
 40 Parham Drive Eastleigh
 SO50 4NU Hampshire
 UNITED KINGDOM
