

# Picus® Electronic Pipette

Catalog No. 106232

## Outstanding Accuracy & Repeatability

- Provides Multi-Dispense Operation with a Single Fill
- Stores Eight User-Defined Pipetting Modes
- Navigates Menu Effortlessly
- Sets Volumes Quickly
- Minimizes Pipetting Errors



The Sartorius **Picus® Electronic Pipette** ensures highly precise and repeatable results. This compact and lightweight pipette has been specially designed to ease the user's workload and to protect the user from repetitive strain injury (RSI).

Perfect  
Companion  
To The  
PAP-8E

## Sophisticated & Ergonomic



Volume Range: 10—300  $\mu$ L

The Multi-Purpose **Adjustment Wheel** Allows Fast One-Handed Operation

The **Electronic Brake** Stops Piston Movement Rapidly and Properly


**Opti-Load** Enables Perfect Tip Sealing for Accurate Delivery from Each Channel

The Intuitive **User Interface** Includes a Choice of Five Language Selections

## Achieve Precise & Consistent Results

### The Multi-Color Backlit Dot-Matrix Display Is Clear & Logically Structured

1. Softkey Functions
2. Current Mode
3. Pipetting Volume & Aliquots in Multi-Dispensing Modes
4. Speed
5. Battery Charging Indicator
6. Arrows Indicating Aspiration & Dispensing



## Additional Pipette Accessories Available

### Picus® Electronic Pipette Tips, Tower Refill



Catalog #106554

- Made with High Quality Polypropylene
- Strict Quality Standards
- Handles 10-350µL of Sample or Reagent

Net Contents: 10 x 96 tips case, 51mm length

### Picus® Electronic Pipette Charging Stand

Catalog #106560

- Charges Battery in Stand
- Prevents Contamination
- Keeps Pipette in Good, Working Order



Dimensions: 380 x 145 x 115 mm

Picus® is a Registered Trademark of Sartorius Lab Instruments GmbH

155 Gibraltar Road  
Horsham, PA 19044 USA

Worldwide: +1 215-441-4000  
USA: 1-800-257-3282  
Fax Worldwide: +1 215-443-8820  
customer.service@biodatacorp.com