



Your PAP-8E Just Got Better!

Upgrade to the Newest, Most User Friendly Software Available for Platelet Aggregation

New Features Include:

- On Demand Optical Calibration Check**
 Optical Calibration may be preformed using purified water. Calibration results printed for regulatory inspections.
- Unlimited User Defined Test Types**
 Define as many tests as needed. Tests may be activated or inactivated according to the user's preference.
- Washed Platelet Sample Type**
 Perform Washed Platelets protocols. Select this test type to auto adjust optical circuit.
- Adjustable Temperature Settings**
 Temperature may be set on user defined tests (32° - 45° C).
- Convenient Database Backup & Recovery**
 Backup or Restore Database to CD Rom or USB Flash Drive.
- Enhanced Password Protection**
 Secure Patient Data through Operator Log-In Time Outs, Periodic Password Expiration and Log-In Failure Lock Out.

For more details, please call 1-800-257-3282 (USA) or 215-441-400 (Worldwide) to speak to our Customer Service Department.



 155 Gibraltar Road
 Horsham, PA 19044 USA

Worldwide: +1 215-441-4000
 USA: 1-800-257-3282
 Fax Worldwide: +1 215-443-8820
 customer.service@biodatacorp.com



An ISO 13485 Registered Company
www.biodatacorp.com



mdi Europa GmbH
 Langenhagener Str. 71
 D-30855 Langenhagen
 GERMANY

Alpha Laboratories
 40 Parham Drive Eastleigh
 SO50 4NU Hampshire
 UNITED KINGDOM

