

The **PAP-8E** Platelet Aggregation Profiler

ROUTINE & SPECIAL AGGREGATION TESTING

RISTOCETIN COFACTOR ACTIVITY ANALYSIS

ANTI-PLATELET DRUG RESPONSES

HEPARIN INDUCED THROMBOCYTOPENIA DETECTION



A COMPLETE TEST SYSTEM

8 Test Channel Aggregometer

All-In-One Wireless Computer

PAP-8E System Software

Electronic Pipette

Microsoft® Office Suite



* Printout Enlarged For Detail

The Gold Standard in LTA Platelet Function Testing

Bio/Data Corporation's PAP-8E is an advanced platelet aggregometer for clinical and research laboratory applications. Exceptional ease of use and standardized test procedures allow for rapid, routine testing of platelet function. The PAP-8E is designed to support detection and diagnosis of hemostasis disorders, including inherited and acquired platelet dysfunctions. The system measures the rate and extent of activation, aggregation, agglutination and inhibition reactions. Comprehensive results include area under the curve, area under the slope, and user defined analyses.

- Eight Independently Operated Test Channels
- Nine Test Results Reported Per Channel
- Standard Micro-Volume (225uL) Samples
- User-Defined Test Templates & Protocols
- Walk-Away Operation
- Touch Screen Control
- On-Demand Optical Calibration Check
- Two Year Limited Warranty

Trusted Accuracy, Best-In-Class Performance, Simple Operation & Rapid Results



