

PICUS[®] 2

Electronic Pipette

Catalog No. 106232

Efficient and Accurate Pipetting for Reliable Performance


Perfect
Companion
To The
PAP-8E

The **PICUS[®] 2** is a single-channel electronic pipette offering volumes ranging from 10 μ L to 300 μ L. Streamline laboratory tasks with its eight readily available application modes. With its ergonomic design, it ensures reliable performance while reducing hand strain, thus ensuring consistent results. Suitable for all skill levels, the **PICUS[®] 2** offers advanced features with a user-friendly interface.

- Effortless Dispensing with Eight Application Modes
- Extensive Storage with Twenty Memory Locations
- Precise Electronic Tip Ejection
- Seamless Bluetooth Connection with Mobile App
- Integrated Service and Calibration Reminder System
- Two-Level Password Protection for Security
- User-Friendly Menu with Color Display
- Improved Battery Life for Extended Use
- Advanced Ergonomic Design



Optimize Your Workflow with the Sartorius Pipetting Mobile App



The Sartorius Pipetting Mobile App is essential for **PICUS[®] 2 Pipette** users. This new App updates the pipette with the latest features, enhances liquid handling precision, and guides workflows step-by-step. Configurable protocols and seamless software integration ensure accuracy and ease of use.


- Connects with Bluetooth or USB
- Allows for Data Transfer and Integration with Laboratory Software
- Enables Remote Monitoring and Control of Pipetting Tasks
- Enhances Workflow Efficiency and Data Management

PICUS® 2 Electronic Pipette

Catalog No. 106232

Advanced Display Layout Provides Real-Time Operating & Program Status

1. Left, Middle, and Right Softkey Functions
2. Current Application Mode
3. Pipetting Volume
4. Aspirating & Dispensing Speed
5. Battery Life Indicator
6. Aspirating or Dispensing Mode Indicator
7. Bluetooth Connectivity



Sophisticated Design for Ergonomic Comfort

The Multi-Purpose **Scrollwheel** Allows Fast One-Handed Operation

The **Electronic Brake** Stops Piston Movement Rapidly and Properly

Opti-Load Feature Enables Perfect Tip Sealing for Accurate Delivery

The Intuitive **User Interface** Includes Language Selections

Additional Picus® 2 Pipette Accessories Available



Picus® & Picus® 2 Pipette Tips, Tower Refill

Catalog # 106554

Dimensions: 51mm length

Net Contents: 1 Case
10 x 96 Tips Per Case




Picus® & Picus® 2 Pipette Charging Stand

Catalog # 106560

Dimensions: 380 x 145 x 115 mm

Net Contents: 1 Stand

Picus® is a Registered Trademark of Sartorius Lab Instruments GmbH

 155 Gibraltar Road
Horsham, PA 19044 USA

Worldwide: +1 215-441-4000
USA: 1-800-257-3282
Fax Worldwide: +1 215-443-8820
customer.service@biodatacorp.com



An ISO 13485 Registered Company

www.biodatacorp.com

REF 106075
106077
106232