

# TRAP-6 amide

Thrombin Receptor Activating Peptide

Catalog No. 106956

**TRAP-6 amide** is a synthetic PAR 1 Antagonist (Thrombin Receptor Activating Peptide). This peptide fragment irreversibly stimulates platelets to *aggregate*, *undergo structural changes*, and to *degranulate*.

## Stimulate Platelet Activation



Net Contents: 5.0mg

Lyophilized for Long Term Stability & Storage

- **INVESTIGATE**  
Platelet Dependent Pro-Coagulant Activity
- **ASSESS**  
Thrombin Initiated Platelet Aggregation
- **CONFIRM**  
Platelet Function & Reactivity In The Presence Of Thienopyridines
- **STUDY**  
Sticky Platelet Syndrome

### SUITABLE FOR:

- \* Platelet Rich Plasma
- \* Flow Cytometry
- \* Gel Electrophoresis
- \* Washed Platelets
- \* Impedance Aggregation
- \* Animal Models

155 Gibraltar Road, Horsham, PA 19044 USA  
**Worldwide:** +1 215-441-4000  
**USA:** 1-800-257-3282  
**Fax Worldwide:** +1 215-443-8820