



AggRecetin®

Ristocetin A-Sulfate

Catalog No. 100970 & 101241

The Gold Standard In Platelet Function Testing

AggRecetin is used in combination with other IVD components, to perform a Routine Platelet Function Test and Ristocetin CoFactor Activity Assay Tests.

AggRecetin induces agglutination of platelets in standardized fixed platelet suspensions and patient plasmas.



AVAILABLE IN TWO SIZES

Catalog No. 100970
Net Contents: 15mg
Includes 2.0mL Diluent

Catalog No. 101241
Net Contents: 100mg
No Diluent Included

Perform RIPA & Ristocetin CoFactor Activity Assays

- Mediates Platelet Adhesion (Agglutination) by Binding to Collagen and Platelets
- Used with Patient Platelet Poor Plasma and Fixed Platelets to Perform the Ristocetin CoFactor Activity Assay
- Causes vWF to Bind to the Gp1B Platelet Receptor
- Causes Ristocetin Induced Platelet Aggregation (RIPA) in Platelet Rich Plasma Samples when added in Two Concentrations
- Determines the Presence and Integrity of the Platelet Glycoprotein (GP) Iba-V-IX Complex and von Willebrand Factor (VWF) Interaction when Performed Using Platelet Rich Plasma



155 Gibraltar Road
Horsham, PA 19044 USA

Worldwide: +1 215-441-4000
USA: 1-800-257-3282
Fax Worldwide: +1 215-443-8820
customer.service@biodatacorp.com

REF
100970
101241



An ISO 13485 Registered Company
www.biodatacorp.com

EC **REP**

mdi Europa GmbH
Langenhagener Str. 71
D-30855 Langenhagen
GERMANY

UK **REP**

Alpha Laboratories
40 Parham Drive Eastleigh
SO50 4NU Hampshire
UNITED KINGDOM

